

Product datasheet for **RR214584**

Dars2 (NM_001034143) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Dars2 (NM_001034143) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Dars2
Synonyms:	RGD1308286
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide
Sequence:

>RR214584 representing NM_001034143
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGTATCTTGGTTCTTGGTTAAATCGTCTAGGCAGGGTTTATCCAGACCCATCGGAAAGACCAAGCAGC
CAATCTGGGGTTCTCTTTCCAGAAGTCTGACGTTGAGTTCACAGAGAGTTCAGAATTCAGTAGCTTTGT
TGCTCGGACCAACACATGTGGAGAACTTCGTTCTCCCACTTAGGCCAAGAAGTCACTCTGTGTGGATGG
ATTCAGTACCGAAGGC AAAACACCTTCTTAGTCTAAGAGATTGCCACGGTCTCGTTCAAATCCTCATT
CCCAGGATGAGTCGGCTGCATCTGTGAGGAGGACTTTGTGTGAAGCCCTGTGGAGTCTGTGGTGAAGT
ATCTGGGACGGTCATCGCCCGACCCCTGGGACAGGAGAATCCGAAAATGCCAACAGGTGAGATTGAAATC
AAAGCTAAACAGCTGAGCTCCTGAATGCCTGCAAAAAGTGCCTTTGAAATTAAGACTTTGTGAAGA
AGACGGAAGCCCTTCGTCTGCAGTATCGGTACTTAGACCTGCGCAGCTCCCAAATGCAGCACAACCTGCG
GCTCAGGTCACAGATGGTCATGAAAATGCGGGAATATCTCTGACTCTGCACGGATTTGTGGATATAGAA
ACCCCAACATTGTTTAAAGAGGACTCCTGGAGGTGCCAAGGAGTTTTTGTAGTACCGTCCAGGGAACCTGGAA
GGTTTTATTCTCTCCCTCAGAGTCTCAACAGTTTAAAGCAGTCTCTGATGGTTGGTGGTTTAGACAGGTA
CTTTCAGGTTGCCCGGTGTATCGAGATGAAGGCTCAAGACCGGATAGGCAGCCAGAGTTTACTCAGATT
GACATAGAGATGTCCTTTGTAGACCAGACTGGCATTGAGCATTAGTTGAGGGTTTGGTGCATTATTCCT
GGCCTGAAGATAAAGACCCTCTAGTGGCTCCTTTCTCCATGACTTTTGTGAGGCATTGGCTACCTA
CGGAAGTATAAGCCAGATACTCGCTTTGGGATGAAGATTGTGGATATCAGTGTGTTCAGAAACACA
GAGATCAGGTTTCTTCAGGATGCACTGGCTAAGCCCAAGGAAGTCAAAGCAATATGTGTTTATGAAG
GAGCAAAGTACTTAAGGAAAGAAGACATTGAGTTCATTTCGAAAGTTTGCAGCTCACCATTTCAGTCAGGA
AGTCTTCTATATTCTGAATGCCAGGAAGAAGTGGAGTTCACCATTTGCTAAATTCATAACGGAGGAG
GAAAGATTGGAGTTAACCAGATTGATGGAGATACAAGAAGACGACATGGTGTGCTGACTGCTGGGCAAC
ATGAGAAAGCATGCTCTTTGTTGGGGAAGTTACGACTAGAGTGTGCAGACCTTTTGGAAACGAGAGGACT
GGCGCTCAGGGACCCAGCATTGTTCTCTTTCTGCGGACTTGGACTTCCCACTCTTCTCGCCAAGGAA
GAGAGTCCACAGAGCTGGAGTCAAGCCACCACCCATTTACTGCTCCTCACCCCGGCGACATTACCTCC
TCTACACTGAGCCGAAAAGGTTCTGGCCAACACTACGACTTGGTTTTAAATGGCAATGAGATAGGAGG
AGGCTCTATCCGATTTCATGATGCACAGCTACAGAGGTACATCCTGGAACATTGCTAAAGGAAGATGTG
AAATTGCTTTCCATCTGCTTCAGGCTTAGACTATGGGGCACCCCCACATGGAGGAATTGCCTTAGGGT
TAGACAGACTGGTGTGCTTGTACAGGGAGCACCCAGCATCAGAGATGTCATAGCCTTCCCTAAATCCTA
CCGAGGACATGACCTTATGAGCAATGCTCCAGACACTGTGTCTCCTGAAGATCTGAAGCCTTACCATATC
CACGTCTCATGGCCAACAGACTCAGAAGAGAGAGCATCCGCCACACCTTCCAAATACTTGAGCTCT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR214584 representing NM_001034143
 Red=Cloning site Green=Tags(s)


MYLGSWLNRLGRGLSRPIGKTKQPIWGSLSRSLTLSSQRVPEFSSFVARTNTCGELRSSHLGQEVTLCGW
 IQYRRQNTFLVLRDCHGLVQILIPQDESAASVRRTLCEAPVESVVRVSGTVIARPLGQENPKMPTGEIEI
 KAKTAELLNACKKLPEIKDFVKKTEALRLQYRYLDRSSQMQHNLRLRSQMVMKMREYLCTLHGFVDIE
 TPTLFRTPGGAKEFLVPSREPGRFYSLPQSPQQFKQLLMVGGLDRYFQVARCVRDEGSRPDRQPEFTQI
 DIEMSFVDQTGIQHLVEGLLHYSWPEDKDPLVAPFPSMTFAEALATYGTDPDTRFGMKIVDISDVFRNT
 EIRFLQDALAKPQGTVKAICVHEGAKYLRKEDIEFIRKFAAHHSQEVLPIFLNARKNWSSPFAKFITEE
 ERLELTRLMEIQEDDMVLLTAGQHEKACSLGKLRLECADLLETRGLALRDPALFSFLWVDFPLFLAKE
 ESPTELESAHHPFTAPHPGDIHLLYTEPEKVRGQHYDLVINGNEIGGGSIRIHDAQLQRYILETLKEDV
 KLLSHLLQALDYGAPPHGGIALGLDRLVCLVTGAPSIRDVIAFPKSYRGHDLMSNAPDTPSPEDLKPYHI
 HVSWPDTSEERASATPSKYLSS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001034143

ORF Size: 1956 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001034143.1](#), [NP_001029315.1](#)

RefSeq Size: 2321 bp
RefSeq ORF: 1959 bp
Locus ID: 304919
UniProt ID: Q3KRD0
Cytogenetics: 13q22
MW: 74 kDa