

Product datasheet for **RR214603**

Ddx49 (NM_001115023) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ddx49 (NM_001115023) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Ddx49
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RR214603 representing NM_001115023
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGC**C

ATGGCGGGGTTTCGCTGAGCTCCGGCTATCGTCGTGGCTGGTGAACAGTGTGACAGCTGGGGTTGAAGC
 AGCCACGCCGGTGCAGCTCGGGTGCATCCCGCCATCCTGGAAGTTCGGGACTGCTTGGTTGTGCAAA
 AACAGGAAGTGGCAAGACAGCAGCCTTTGACTTCCCATCCTGCAGAAGTGTCTGAGGACCCTTATGGC
 ATCTTCTGCCTGGTTTTGACACCACTAGAGAGCTGGCCTACCAGATAGCTGAGCAGTTCGAGTGCTGG
 GGAAGCCTCTGGGCCTAAAGGACTGCATCATCGTTGGTGGCATGGACATGGTGACGCAGGCCCTGGAGCT
 CTCTCGGAAGCCACAGTGGTTATTGCCACCCCTGGGCGCCTGGTGACCACTGCGCAGTTCACCACT
 TTCAACATGAAGAAAATCAGATTCTTGGTGTGGACGAGGCCGATAGGCTGCTGGAGCAGGGTGCACCG
 ACTTCACTGCAGACCTGGAGACCATCCTGTGAGCTGTTCCAGCACGCAGGCAACCCTCCTGTTCACTGC
 CACACTGACGGACACACTCAAGGAGCTGCAGGGCCTGGCTACCAACCAGCCCTTCTTTGGGAGGCACAG
 GCCACGGTTCGCATGGTGGAGCAGCTTGACCAGCGCTATCTGCTGGTGCCTGAGAAGGTGAAGGACGCCCT
 ACCTGGTTCACCTGGTCCAGACCTTCCAGGACCAGCTGGAGGACTGCTCCATCATCATCTTCACCAATAC
 CTGCAAGACCTGCCAGATCCTTTGTATGATGCTTCGAAAGTTCAACTTCCCCACAGTGGCACTGCACTCC
 ATGATGAAGCAGAAAGAGCGCTTTGAGCCCTGGCCAGATTCAAGTCCAGCACCTATCGGATTCTTATTG
 CTACTGACGTGGCCTCCAGAGCCTGGACATCCCCACTGTGCAGGTGGTCATCAACCACAACACCCTGG
 GCTCCCCAAGATCTACATCCATAGGGTCGGTCGCACAGCCCGTGCAGGGAGGCAGGGACAGGCCATCACA
 CTGGTGACGCAGTACGACATCCACCTGTTGCATGCCATCGAGGAGGAGATCAAACAGCAGCTGGCTGAAT
 TGGCAGTGGAGGAGGCCAGGTCTGCAGATCCTCACACAGGTCAACGTGGTGGCAGGGAGTGTGAAAT
 CAAACTGGAGGCTTCGCATTTTGTGATAAAGAGGAGATCAACAAGCGGAAACAGATGATCCTGGAGGGG
 AAGGACCCGACCTGGAGGCCAAGCGCAAAGCTGAGCTGGCCAAGATTAAGCAACAGAACCAGCGGCTTCA
 AGGAGAAAGTGGAGCAGACGCTTCGGCAGCAGAAAGCAGGGAGCACAGGCCGGCGGTGCGTCTCCAG
 ACCTCGTCCCCAAGAGCCAGCCAGGGCTGAGGCTCAGAAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RR214603 representing NM_001115023
 Red=Cloning site Green=Tags(s)


MAGFAELRLSSWLVEQCRQLGLKQPTPVQLGCIPAILEGRDCLGCAKTGSGKTAAFVLPILQKLSDEPYG
 IFCLVLTPTRELAYQIAEQFRVLGKPLGLKDCIIVGGMDMVTQALELSRKPHVVIATPGLADHLRSSNT
 FNMKKIRFLVLDEADRLLEQGCTDFTADLETILSAVPARRQTLTFSATLTDTLKELQGLATNQPFWEAQ
 ATVRMVEQLDQRYLLVPEKVKDAYLVHLVQTFQDQLEDCSIIFNTCKTCQILCMMLRKFNFTVALHS
 MMKQKERFAALARFKSSTYRILIATDVASRGDIPTVQVVINHTPGLPKIYIHRVGR TARAGRQQAIT
 LVTQYDIHLLHAIEEIIKQQLAELAVEEAQVLQILTQVNVVRRECEIKLEASHFDEKKEINKRKQMILEG
 KDPDLEAKRKAELAKIKQNRFRKEKVEQTLRQQKAGSTGRRSRPPRPQEPARAEAQN

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Restriction Sites:


Sgfl-MluI

Cloning Scheme:



ACCN:	NM_001115023
ORF Size:	1440 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001115023.1 , NP_001108495.1
RefSeq Size:	1555 bp
RefSeq ORF:	1443 bp
Locus ID:	290660
Cytogenetics:	16p14
MW:	54.3 kDa

Product images:



Circular map for RR214603