

## Product datasheet for **RR214646L4V**

### Fahd2a (NM\_001134834) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Fahd2a (NM_001134834) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Fahd2a
Synonyms:	Fahd2; RGD1563674
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001134834
ORF Size:	939 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR214646).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001134834.1</a> , <a href="#">NP_001128306.1</a>
RefSeq Size:	1257 bp
RefSeq ORF:	942 bp
Locus ID:	296131
UniProt ID:	<a href="#">B2RYW9</a>
Cytogenetics:	3q36
Gene Summary:	May have hydrolase activity.[UniProtKB/Swiss-Prot Function]



[View online »](#)