

## Product datasheet for **RR214647L3V**

### Pax6 (NM\_013001) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Pax6 (NM_013001) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Pax6
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_013001
ORF Size:	1266 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR214647).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_013001.2</a> , <a href="#">NP_037133.1</a>
RefSeq Size:	2191 bp
RefSeq ORF:	1269 bp
Locus ID:	25509
UniProt ID:	<a href="#">P63016</a>
Cytogenetics:	3q33
Gene Summary:	homeodomain (HD)-containing transcription factor; may play a role in the differentiation of the hindbrain somatic (SM) motoneurons and V1 interneurons in the hindbrain and/or spinal cord [RGD, Feb 2006]



[View online »](#)