

Product datasheet for **RR214661**

Loxl1 (NM_001012125) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Loxl1 (NM_001012125) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Loxl1
Synonyms:	MGC105563
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>RR214661 representing NM_001012125
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGCTCTGGCCGGAGCCGGCAGCCAGCTGGGGACCTGGTGTGGGGAGCCTGCCTGTGCGTCTCGTGC
 ACGGGCAGCAAGAGCAGCCAGGACAGGGCTCGGACCCAGGACGCTGGCGGCAGCTGATCCAGTGGGAGAA
 CAACGGGCAGGTGTACAGCTTGTCAATTCGGGCTCTGAGTACGTGCCTGCGGGGCCCAAGAGGGGAG
 AGTAGCTCCAGGGTGTACTGGCAGGCGCGCCCAAACGTGCGAGCGGCGCAGTCAGGGAGGACCCCGGC
 GTCGGCAGGCGCCGAACCTGCCCTGCCCGGACGCGTGGGTTCGGACTGTGCGCGGCCAAACGCGGCA
 CCCGTTTCGGCTTCGGCCAGGTGCCAGACAACCTGGAGAGAGGTGGCCGTGCGGGACAGTAACGGGATGGCT
 CGGGCCCGCACCTCCGTCTCCAGCAACGGCACGGGGCTCCGCTCCTCCTCGGTTTCGGCTCGGCT
 TCGCCACCACTTATCGCCAGCCATCCTCTACCCGAGCAGTTCCTTATCCGAGGCGCCCTTCGTGAA
 CCAATATGAGAACTATGACCCCGCTGCGGACCTATGAACAGGGCTACGTGTACTACCGGGTGGGGC
 GGAGGCATGGGCGCAGGGGCGGCGGCGGACGCTGAGCCTCGGCTGGGGTCATCTATCCCTCCAGCCCC
 GGGCACGCTACGAGGAGTACGGCGGTGGGGGCGAGGAACAGCCCGAGTACCCGGCGCAGGGCTTCTACCC
 AGCCCCAGAGAGGCCCTATGCGCCCCAGCCCCAGCCTCAGCCCCAACCTCAGCCGACGCTCAGCCCCAA
 CCTTCCGACGGCTTGACCGTCTTACTCGCATAGCCTGTACAATGAAGGCACCCCGGGATTTCGAGCAGG
 CCTACCTGATCCAGCACCGACGTGTCGAGATCCCCGCCGGGGCCGGTGGCACCTACGGTGGCCGTGG
 CGATCCCCGTCTCGGCTGGTACCCGCCCTATGCAGCAATGTGCCTCCGAGGCCATGTCCCTCCGCGT
 GCGGTGGAGCCGACGCCCATTCGCGTGTGGAACACCTTACTTGCCTGTGCGAAACTCTGATCGGC
 CCTCCAGGGTGGGAGAGGAACGGTGCCAGCAGGGCCGTCTCAGCGTGGGTAGTGTGTACCGACCTAA
 CCAGAATGGTTCGCGGCTCCCTGACTTAGTCCAGATCCAACTACGTGCAGGCATCCACATATGTGCAG
 AGAGCCCATCTGTATTCCTTGCCTGTGCGGCAGAGGAGAAGTGCCTGGCCAGCACAGCCTATGCTCCTG
 AGGCCACCGACTATGACCTCCGAGTGTGTACGATTCCCCAGCGTGTGAAAAACAGGGCACGGCGGA
 CTTCTCCCAACCGGCCACGCCACCTGGGAGTGGCAGTGTGACCAACATTACCACAGCATGGAC
 GAGTTCAGCCACTATGACCTACTTGTGATGCCACGGGCAAGAAAGTTGCCGAGGGACACAAGGCCAGCT
 TTTGCCTGGAGGACAGCACCTGTGACTTCGGCAACCTCAAGCGCTATGCGTGCACCTCACACACGAGGG
 CCTGAGTCCAGGCTGCTATGACACCTACAATGCGGACATCGACTGCCAGTGGATGCATAACAGATGTG
 CAGCCTGGAACTACATCCTCAAGGTGCAGTGAACCCCAAGTACATTGTTCTGGAGTCTGACTTACCA
 ACAACGTGGTGTAGATGCAACATCCACTACACAGTTCGCTCTCTACAACAACTGCAAAATCGTCCA
 GTCC

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
 TGGATTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:


>RR214661 representing NM_001012125
 Red=Cloning site Green=Tags(s)

MALAGAGSQLGTLVWGAACLVLVHGQEQPGQSDPGRWRQLIQWENNGQVYSLNLSGSEYVPAGPQRGE
 SSSRVLLAGAPQTSQRRSQGGPRRRQAPNLPLPGRVGS DTVRQTRHPFGFQVDPDNWREAVGDSNGMA
 RARTSVSQQRHGGSSASSVSASAFATTYRQPSSYPQQFPYPQAPFVNQYENYDPASRTYEQGVVYRAGG
 GGMGAGAAAAAVASAGVIYPFQPRARYEYGGGEEQPEYPAQGFYPAPERPYAPQPQPQPQPQPQPQ
 PSDGLDRRYSHSLYNEGTPGFEQAYPDPSTDVSIIPAGAGGTGAAAGDPRLGWYPPYAANVPPEAYVPPR
 AVEPQPPFRVLEPPYL PVRNSDAPSQGGERNGAQQGRLSVGSVYRPNQNGRGLPDLVPDPNYVQASTYVQ
 RAHL YSLRCAAEKCLASTAYAPEATDYDLRVLLRFPQRVKNGTADFLPNRPRHTWEWHSCHQHYHSM
 EF SHYDLLDASTGKKVAEGHKASF CLEDSTCDFGNL KRYACTSHTQGLSPGCYDTYNADIDCQWIDITDV
 QPGNYILKVHVNPKYIVLESDF TNNVVR CNIH YTRFVSTTNCKIVQS


SGP**TRTRPLEQKLI**SEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-RsrII

Cloning Scheme:



Plasmid Map:



ACCN: NM_001012125

ORF Size: 1824 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001012125.1 , NP_001012125.1
RefSeq Size:	2491 bp
RefSeq ORF:	1827 bp
Locus ID:	315714
Cytogenetics:	8q24
MW:	66.6 kDa