

Product datasheet for **RR214662L4V**

Kcnj5 (NM_017297) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Kcnj5 (NM_017297) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Kcnj5
Synonyms:	MGC93525
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_017297
ORF Size:	1257 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR214662).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_017297.2 , NP_058993.1
RefSeq Size:	3335 bp
RefSeq ORF:	1260 bp
Locus ID:	29713
UniProt ID:	P48548
Cytogenetics:	8q21
Gene Summary:	inwardly-rectifying potassium channel which is controlled by G proteins [RGD, Feb 2006]



[View online »](#)