

## Product datasheet for **RR214671**

### **Ptk6 (NM\_001108968) Rat Tagged ORF Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Ptk6 (NM\_001108968) Rat Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Ptk6  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RR214671 representing NM\_001108968  
**Red=Cloning site Blue=ORF Green=Tags(s)**

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGTGTCTTTGGACAAGGCTCACCTGGGTCCTAAGTATGTGGCCCTCTGGGACTCAAAGCACGAACGG  
ACGAGGAGCTGAGCTTTACGGCAGGAGACCTCCTCCATGTTACCAGGAAGGAGGAACAGTGGTGGGC  
CACCTGCTGGATGCGGAAGGCAAGGCCTTGGCTGAGGGCTACGTGCCCCACAACCTACCTGGCTGAGAAG  
GAAACTGTGGAGTCTGAACCGTGGTTCTTTGGTTGCATCTCCCGCTCAGAGGCCATGCACAGGCTGCAAGT  
CTGAGGACTACCCGAAGGTACCTTCTGATCAGAGTCAGCCAGAAGCCAGGAGCAGACTATGTCCTCTC  
CGTGTGGGATGCTGAGGCTGTTGCACATTACAGGATCTGGAGAACAGTGCAGGAGGCTGCACCTGAAT  
GAGGTGGTATCCTTCTCCAGTCTGTCTGAGCTTGTGAACACCACAAGACCCACAACCTGTCCCCTGGCC  
TACAGCTCTCCATGGCCTGCTGGAAGCTCAAACTGAGCCCTGCCCCACTGGGATGACTGGGAGAGGCC  
GAGGGAGGAGTTCACGCTCTGTAAGAAGCTGGGGTCTGGCTACTTTGGGAGGTCTTTGAAGGGCTCTGG  
AAAGCCCTGGTCCGTGTGGCTGTGAAGGTGATCTCTAGAGACAATCTCCTGCACCAGCACACCTTCCAGG  
CTGAGATTCAGCCATGAAGAAGCTGCGACACAAGCACATCCTGTCGCTGTACGCTGTGGCAACTGCAGG  
GGACCCAGTCTATATCATCACGGAGCTCATGCCAAGGGGAGCCTGCTGCAGCTATTGCGGGACTCGGAT  
GAGAAAGACCTGCCTATTTTGGAACTGGTGGACTTTGCATCACAGTTGCCGAGGGCATGTGCTACCTGG  
AATCTCAGAATTACATTCACCGTACCTGGCTGCAAGGAACGTTCTCGTTACTGAGAACAATCTCTGCAA  
AGTGGGGGACTTTGGGCTGGCCAGGCTCGTCAAGGAGGACATCTACCTTTCCCGTGACCACAACGTCCTCC  
TACAAATGGACAGCACCTGAGGCACTTTCCCGAGGGCATTACTCCATCAAATCGGATGTCTGGTCTTTTG  
GTGTTCTCCTTCAATAATTTTTCAGCAGGGGCGAGATGCCCTACCCAGGCATGTCCAATCACGAAACCTT  
CCTGAGGGTGGATGCCGGCTACCGCATGCCCTGCCCTTGGAGTGCCACCCAACATACACAAGCTGATG  
CTTAGCTGTTGGAGCAGAGACCCTAAGCAGAGACCCTGCTTCAAAGACCTGTGTGAGAACTCTCAGGTA  
TCACTAGGTATGAGAACCTGGTC

**ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT**  
**ACAAGGATGACGACGATAAGGTTTAA**




**Protein Sequence:** >RR214671 representing NM\_001108968  
 Red=Cloning site Green=Tags(s)

MVSLDKAHLGPKYVGLWDFKARTDEELSFQAGDLLHVTRKEEQWWATLLDAEGKALAEGYVPHNYLAEK  
 ETVESEPWFFGCISRSEAMHRLQSEDPKGTFLIRVSQKPGADYVLSVWDAEAVRHYRIWRNSAGRLHLN  
 EVVSFSSLSELVNYHKTHNLSPGLQLSMACWKLKTEPLPHWDDWERPREEFTLCKKLGSGYFGEVFEGLW  
 KGLVRVAVKVISRDNLLHQHTFQAEIQAMKKLRHKHILSLYAVATAGDPVYIITELMPKGSLLQLLRDSD  
 EKDLPILELVDFASQVAEGMCYLESQNYIHRDLAARNVLVTENNLCKVGFGLARLVKEDIYLSRDHNVP  
 YKWTAPEALSRGHYSIKSDVWSFGVLLHEIFSRGQMPYPGMSNHETFRLVDAGYRMPCLPCPPNIHKLM  
 LSCWSRDPKQRPCKDLCEKLSGITRYENLV


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**Plasmid Map:**



<b>ACCN:</b>	NM_001108968
<b>ORF Size:</b>	1353 bp
<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001108968.1</a> , <a href="#">NP_001102438.1</a>
<b>RefSeq Size:</b>	2142 bp
<b>RefSeq ORF:</b>	1356 bp
<b>Locus ID:</b>	366275
<b>Cytogenetics:</b>	3q43
<b>MW:</b>	52 kDa