

Product datasheet for **RR214685**

Ccr2 (NM_021866) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Ccr2 (NM_021866) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Ccr2
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR214685 representing NM_021866
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGAAGACAGTAATATGTTACCTCAGTTCATCCATGGCATACTATCAACATCTCATTCTCTATTCCAA
GAAGTATCCAAGAGCTTGATGAGGGGGCCACCACCCGATGACTATGATGATGGTGAACCTTGTCCATAA
AACCAGTGTGAAGCAAATTGGAGCTTGGATCCTGCCCCACTCTACTCCCTGGTATTCATCTTTGGTTTT
GTGGGCAACATGTTGGTCATTATAATTCTGATAAGCTGTAAAAAGCTGAAGAGCATGACTGATATCTACC
TGTTCAACCTGGCCATCTCTGACCTGCTCTTCTGCTCACACTCCCATTCTGGGCTCACTATGCTGCAAA
TGAGTGGGTCTTTGGGAATATAATGTGCAAATTATCACAGGGCTTTATCACATTGGGTATTTGGTGGAA
ATCTTCTCATTATCCTCCTGACAATTGATAGATATTTGGCTATTGTCCATGCTGTCTTTGCTTTAAAAG
CCAGGACAGTTACCTTTGGGTAATAACAAGTGTAGTCACTTGGGTGGTGGCTGTGTTGCCTCTCTACC
AGGAATCATATTTACTAAATCTGAACAAGAAGATGATCAGCATACTTGTGGCCCTATTTTCCAACAATC
TGGAAGAATTTCCAAACAATAATGAGGAATATCTTGAGTTTGATCCTGCCCTACTTGTGCATGGTCATCT
GCTACTCAGGAATCCTCCACACCCTGTTTCGCTGTAGGAATGAGAAAAAGAGGCATAGGGCTGTGAGGCT
CATCTTTGCCATCATGATTGTCTACTTTCTCTTCTGGACTCCATACAATATTGTTCTCTTCTGACCACC
TTCCAGGAATCTTGGGAATGAGTAAGTGTGTTGACATGCACTTAGACCAGGCCATGCAAGTGACAG
AGACTCTTGGAAATGACACACTGCTGCGTTAATCCTATCATTTATGCCTTTGTTGGTGAAGAAGTCCGAAG
GTATCTCTCCATATTTTTCAGAAAGCACATTGCCAAAAATCTCTGCAAACAATGCCCAAGTTTCTATAGG
GAGACAGCAGACCAGTGAGCTCAACATTTACCCTTCTACTGGGAGCAAGAAGTCTCAGTTGGGTTG

ACGCGTACGCGGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR214685 representing NM_021866
 Red=Cloning site Green=Tags(s)

MEDSNMLPQFIHGILSTSHSLFPPRSIQELDEGATTPYDYDDGEPCHKTSVKQIGAWILPPLYSLVFIFGF
 VGNMLVIIILISCKKLKSMTDIYLFNLAI SDLLFLLTLPFWAHYAANWVFGNIMCKLFTGLYHIGYFGG
 IFFIILLTIDRYLAIVHAVFALKARTVTFGVITSVVTWVAVFASLPGIIFTKSEQEDDQHTCGPYFPTI
 WKNFQTIMRNILSLILPLLVMVICYSGILHTLFCRNEKKRHRAVRLIFAIMIVYFLFWTPYNIIVLFLTT
 FQEF LGMSNCVDMHLDQAMQVTE TLGMTHCCVNP IYAFVGEKFRRLYSIFFRKHI AKNLCKQCPVFYR
 ETADRV SSTFTPSTGEQEVSVGL


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_021866

ORF Size:	1119 bp
OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_021866.1 , NP_068638.1
RefSeq Size:	1122 bp
RefSeq ORF:	1122 bp
Locus ID:	60463
UniProt ID:	O55193
Cytogenetics:	8q32
MW:	42.8 kDa
Gene Summary:	mediates leukocyte migration and activation [RGD, Feb 2006]