

Product datasheet for **RR214785L3V**

Arhgap27 (NM_198759) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Arhgap27 (NM_198759) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Arhgap27
Synonyms:	Arhgap15; CAMGAP1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_198759
ORF Size:	2607 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR214785).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_198759.1 , NP_942054.1
RefSeq Size:	2610 bp
RefSeq ORF:	2610 bp
Locus ID:	303583
UniProt ID:	Q6TLK4
Cytogenetics:	10q32.1
Gene Summary:	has GTPase activator activity toward Rac1 and Cdc42; plays a role in receptor mediated endocytosis [RGD, Feb 2006]



[View online »](#)