

## Product datasheet for **RR214801L3V**

### Slc25a4 (NM\_053515) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Slc25a4 (NM_053515) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Slc25a4
Synonyms:	Ant1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_053515
ORF Size:	894 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR214801).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_053515.1</a> , <a href="#">NP_445967.1</a>
RefSeq Size:	1289 bp
RefSeq ORF:	897 bp
Locus ID:	85333
UniProt ID:	<a href="#">Q05962</a>
Cytogenetics:	16q11
Gene Summary:	catalyzes the exchange of ADP and ATP across the mitochondrial inner membrane [RGD, Feb 2006]



[View online »](#)