

## Product datasheet for **RR214811L4V**

### Maged1 (NM\_053409) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Maged1 (NM_053409) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Maged1
Synonyms:	Nrage
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_053409
ORF Size:	2325 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR214811).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_053409.2</a> , <a href="#">NP_445861.1</a>
RefSeq Size:	2695 bp
RefSeq ORF:	2328 bp
Locus ID:	84469
UniProt ID:	<a href="#">Q9ES73</a>
Cytogenetics:	Xq22
Gene Summary:	interacts with the p75 neurotrophin receptor; involved in mediating nerve growth factor-dependent apoptosis via a JNK-dependent pathway [RGD, Feb 2006]



[View online »](#)