

Product datasheet for RR214818L3

Gys1 (NM_001109615) Rat Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Gys1 (NM_001109615) Rat Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	Gys1
Synonyms:	MGC156743
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR214818).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.



[View online »](#)

Locus ID: 690987
UniProt ID: [A2RRU1](#)
Cytogenetics: 1q22
Gene Summary: Transfers the glycosyl residue from UDP-Glc to the non-reducing end of alpha-1,4-glucan.
[UniProtKB/Swiss-Prot Function]