

Product datasheet for **RR214820**

Myoz2 (NM_001106469) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Myoz2 (NM_001106469) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Myoz2
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RR214820 representing NM_001106469
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGACACGCACAATCTTGCCATGCAGAATGGGAGAATTGACGGAAGTAACCTGGAAGGTGGTTCACAGC
 AAGCTCCCTCAACTCCTCCTAACACCCCGGATCCACGAAGCCCCCAATCCGGAGAACATTGCACCAGG
 ATATCTGGACCCTGAAGGAAATTCCTCCTGAAAGATTAACACGACGGCTGTTCCTAAGTACTACCAG
 TCTCCCTGGGAGCAGGCCATTGGCGACGATCCAGAGCTCCTGGAAGCTTTGTACCCGAAACTTTTCAAGC
 CTGAAGGAAAGGCAGAACTGCCAGATTACAGGAGCTTTAACAGAGTTGCCACTCCATTTGGAGGTTTTGA
 GAAAGCATCAAAAATGGTTAAATTCAAGGTTCCAGACTTCGAACGACTGCTGCTAACAGATCCCAGGTTT
 TTGGCCTTTGCCAATCCTCTTTTCGGCAGACGATGCTTTAATAGGACACCAAAGGGTGGATATCTGAGA
 ATATTCCTGTCGTGATCACAACCTGAGCCTGCAGAAGATGTCACTATGCCAGAGTCAGATGACTTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR214820 representing NM_001106469
 Red=Cloning site Green=Tags(s)

MDTHNLAMQNGRIDGSNLEGGSQAPSTPPNTPDPRSPNPENIAPGYSGLKEIPPERFNNTAVPKYYQ
 SPWEQAIGDDPELLEALYPKLFKPEGKAEIPDYRSFNRVATPFGGFEKASKMKVFKVPDFERLLLTDPRF
 LAFANPLSGRRCFNRTPKGWISENIPVVITTEPAEDVTMPESDDL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_001106469

ORF Size: 555 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001106469.1](#), [NP_001099939.1](#)

RefSeq Size: 951 bp


RefSeq ORF: 558 bp

Locus ID: 295426

Cytogenetics: 2q42

MW: 20.6 kDa

Product images:



Circular map for RR214820