

Product datasheet for **RR214897L3V**

Sh3kbp1 (NM_053360) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Sh3kbp1 (NM_053360) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Sh3kbp1
Synonyms:	CIN85; Seta
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_053360
ORF Size:	1995 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR214897).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_053360.2 , NP_445812.1
RefSeq Size:	3310 bp
RefSeq ORF:	1998 bp
Locus ID:	84357
UniProt ID:	Q925Q9
Cytogenetics:	Xq14
Gene Summary:	binds apoptosis-linked gene 2 (ALG-2) interacting protein 1; involved in regulation of apoptosis in astrocytes [RGD, Feb 2006]



[View online »](#)