

Product datasheet for **RR214905L3V**

Rhd (NM_022505) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Rhd (NM_022505) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Rhd
Synonyms:	Rh; Rhcd
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_022505
ORF Size:	1266 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR214905).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_022505.2 , NP_071950.1
RefSeq Size:	1509 bp
RefSeq ORF:	1269 bp
Locus ID:	60414
UniProt ID:	O88298
Cytogenetics:	5q36
Gene Summary:	blood group gene; produces antigens on the surface of erythrocytes [RGD, Feb 2006]



[View online »](#)