

Product datasheet for **RR214906L4V**

Chrm4 (NM_031547) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Chrm4 (NM_031547) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Chrm4
Synonyms:	M4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_031547
ORF Size:	1434 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR214906).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_031547.1 , NP_113735.1
RefSeq Size:	2768 bp
RefSeq ORF:	1437 bp
Locus ID:	25111
UniProt ID:	P08485
Cytogenetics:	3q24
Gene Summary:	a G-protein coupled receptor for acetylcholine; important for mediating the effects of acetylcholine in the nervous system, exocrine glands and smooth muscle [RGD, Feb 2006]



[View online »](#)