

Product datasheet for **RR214942**

Tkfc (NM_001039031) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Tkfc (NM_001039031) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Tkfc
Synonyms:	Dak; RGD1311026
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>RR214942 representing NM_001039031
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGAGCTCCAAGAAGATGGTGAATCTGTGGAAGGGTGTGCTGGTATGCCCTTGTGGGTTTGTTCCT
 GTAACCTGACTTGCAGCTCCTGCAAGGGTACCGTGTAGCCCTACGTTCTGACCTGGACAGCCTCAAGG
 CCGGGTGGCACTTCTGTAGGTGGAGGCTCAGGCCATGAGCCTGCCATGCCGGTTTTATCGGAAAGGA
 ATGCTGACAGGTGCATTGCAGGAGCAGTGTTCCTCTCCTGTGTGGGACAGCATCTGGCAGCCATTA
 GAGCTGTGGCCAGGCAGGCACAGCCGGCACCTCCTCATCGTGAAGAACTACACTGGGGATCGGCTCAA
 CTTTGGACTTGCCATGGAGCAGGCCAAGGCCGAGGGCATCTCTGTGGAGATGGTGGTAATTGAGGATGAC
 AGCGCCTTCACTGTCCTGAAGAAGGCAGGCCGGCTGGCCTGTGTGCACAATACTTATCCACAAGGTAG
 CAGGTGCTCTAGCTGAGGAAGGTATGGGCTGGAGGAGATCACAAGAAGGTGAGCGTGATGCCAAGGC
 CATAGGTAATCTGGGGTGGAGCTGTCTCCCTGCAGTGTCCCTGGTACCAAACCCACCTTTGAGCTTGCA
 GCTGATGAAATGGAGCTAGGCCCTGGGGATCCATGGGGAAGCCGGTGTTCGTCGGATAAAGCTGGTACCTG
 TTGATCAGATTGTCACTTTATGCTTGACCACATGACAGACACCTCCAATATATCCCATGTGCCTGTGAA
 GTCAGGCTCTTCACTGGTGTGATGGTCAACAATCTGGGTGGCCTGTCTTTCCTGGAAGTGGGCATCATA
 GCTGATGCTGCCATACGCCTGCTGGAGGGCCGTGGGGTGAAGGTCGCCCTGCCCTGGTGGGTACCTTCA
 TGTGACACTAGAGATGCGTGGTGTTCCTTACTTTGATGCTTGTGGATGAACCCCTGCTGAAGCTGAT
 AGATGCTGAACTAATGCAAAAGCCTGGCCTCACATGTCCAAGGTCTCTGTGACTGGAAGGAATCGGATT
 CGGGCAGCTCCCACAGAGCCTGCAGAAGCCCTGAAGCCACTGCAGCAGGAGGTGTAGCATCAAAGCAGA
 TGACACTTGTGCTGGATCGGATATCCACCACGCTCATCGGACTGGAGGAGCACCTGAATGCCCTGGACAG
 GGCTGCTGGTGTGGGGACTGTGGTTCTACCCATAGCCGTGCTGCCAAAGCCATTACAGGGTGGTTAAAG
 GAAGGCCCAACTCCAGCCAGCCCTGCCAGGTGCTCTCAAACCTGTCTGTCTGCTGGAGAAGATGG
 GAGGCTCATCTGGGGCGCTCTACGGCCTGTTCTGACTGCAGCTGCCAGCCTCTCAAGGCCAACACTGA
 CCTCCCAGCCTGGTCTGCTGCTATGGATGCCGGCCTAAAGGCCATGCAGAAGTATGGAAAGGCTGCACCA
 GGAGACAGGACAATGCTGGATTCTCTGTGGCAGCAGCACAGGAGCTCCAAGCCTGGAAGAGCCAGGGG
 CCAGTCTCCTCCCAGTCTGACTAAAGCAGTCAAAGTGCCGAAGCTGCAGCCGAGGCCACCAAGAACAT
 GGAAGCTGGTGTGGGAGAGCGAGTTACATCAGCTCTGCACAACCTAGATCAGCCAGACCCTGGAGCAGTT
 GCGGCTGCAGCCATCTCCGTGCCATCCTGGAGGTCTTGCAGACGAAGGCAGT

ACGCGTACGCGGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RR214942 representing NM_001039031
 Red=Cloning site Green=Tags(s)


MSSKMKVNSVEGCAGDALAGFVACNPDLQLLQGYRVALRSDLSLKGRVALLSGGSGHEPAHAGFIGKG
 MLTGVIAGAVFASPAVGSILAAIRAVAQAGTAGTLLIVKNYTGDRLNFLAMEQAKAEGISVEMVIEDD
 SAFTVLKAGRRRLCGTILIHKAVAGALAEEMGLEEITKKVSVIAKAIGTLGVSLSPCSVPGTKPTFELA
 ADEMELGLGIHGEAGVRRIKLVPVDQIVTLMLDHMTDTSNISHVPVKS GSSVLMVNNLGGLSFLELGI
 ADAAIRLLEGRGVKVARALVGTFM SALEMRGVSLTLMLVDEPLLKLIDAETNAKAWPHMSKVSVTGRNRI
 RAAPTEPAEPEATAAGGVASKQMTLVLDRISTTLIGLEEHLNALDRAAGDGDGCGSTHSRAAKAIQGWL
 EGPTPASPAQVLSKLSVLLLEKMGSSGALYGLFLTAAQPLKANTDLPAWSAAMDAGLKAMQKYGKAAP
 GDRTMLDSLWAAQELQAWKSPGASLLPVLTKAVKSAEAAAEATKNMEAGAGRAS YISSAQLDQDPDPAV
 AAAAIFRAILEVLQTKAA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-Mlul

Cloning Scheme:



Plasmid Map:



ACCN: NM_001039031
 ORF Size: 1734 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001039031.1 , NP_001034120.1
RefSeq Size:	2403 bp
RefSeq ORF:	1737 bp
Locus ID:	361730
UniProt ID:	Q4KLZ6
Cytogenetics:	1q43
MW:	59.4 kDa
Gene Summary:	Catalyzes both the phosphorylation of dihydroxyacetone and of glyceraldehyde, and the splitting of ribonucleoside diphosphate-X compounds among which FAD is the best substrate (PubMed:16289032). Represses IFIH1-mediated cellular antiviral response (By similarity). [UniProtKB/Swiss-Prot Function]