

Product datasheet for **RR214955**

Slc44a3 (NM_001013914) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Slc44a3 (NM_001013914) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Slc44a3
Synonyms:	RGD1305808
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RR214955 representing NM_001013914
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGGTTACTCGTGGTGGCCGGGGCGGCAGGAAGACTCCTCTTTGGATATGACAGCTTTGGCAATGTGT
 GTGGCAAGAGGAACTCTCCGGTGAAGGCGCCCCCTCTCTGGCAGGACATGACTCTAAAAAACATGT
 GTTCTTTATGAATGCCTGCAATCTGGAAGTCAAAGACAGGGGACGTGGTCCCATGGCCCTCTGTGTATCC
 AGCTGTCCCGAGAAGCAGCTGGACACCCTGGAAGAGGTCCAGCTTTTCGCCAACATCAATGGCTCTTTCC
 TGTGTGTATAGTCTGAATTCCTCAACTACACTCAGAGTCCAAGTGCAGATAGACTGTGTCCCAGGCT
 CCCGGTTCCTCCGAGCAAGCCGTTCCCCTATTTAACCGGTGCATCCCTCAGACCCCGAGTGCTATTCC
 TTGTTTGCATCGGTCTTGCTGAACGATGCCGACGCTCCATCGGATTCTAAGTGAATCATGGCAGGAA
 GAGACACCATCCTCGGCTTGTGTGCTTTGTATTCGCCTTGTCTCTGGCCATGTTGTTACCTCCGATT
 CATCTCAACGCTCCTGGTCCACATTATCATCTCCTTGGTCATTCTGGGATTGCTGTTTGTCTGTGGTGTG
 TTCTGGTGGCTGTATTATGACTACACCAATGACCTCAGTACAGAATTGGACACAGAAAAGGAGAATATGA
 AGTGCATGCTGGCCTTGTCTCATCACCACCGTCGTGACAGTCTTCTCCTGGCCCTGATTTTACCCT
 CAGAAAGAGAGTAAAATTGACGGTTGAGCTTCTCCGAGTCACAAATAAAGCAATCAGCAGATGTCCTTTC
 CTGCTGTTGCAGCCACTGTGGACATTTGCCATCCTCGTTTTCTCTGGGTCTCTGGGTGGCTGTGCTGC
 TGAGTCTGGGAACCGCAGGACTGCCAGGTGATGGAAGTGGCCAAGTGAATATAAACCTCTCTCAGG
 CATTCCGTATCTGTGGTGGTACCATTTAATTGGCTCATCTGGACTAGTGAATTCATCCTTACGTGCCAG
 CGGATGACTGTAGCAGGAGCGATGGTCACTTGTATTTCAACAGGAATCAAAATGACCCTCCTGCCGTC
 CAATCCTTTTATCTCTCCACACTTTTCTGCTACCATCAAGAACTGCAGTGAAGGGTCACTTTTACT
 CACTGTGACAAGGATCCCAAGAGTCATTTTCATGTACATTTATAGCACCGTGAAAGAGCGGCACAGTGCT
 TGGCCAAGAGGTGAATTCAGATGCAGCTACTGTGGCTTTGGTGTCTACCAATTATCCATACCATCTCA
 ATCAGGACCGTATGCCGACGGCCATTAATGGGACAGACTTTTGTACATCAGCAAAAGACGCGCACAC
 AATCATATCCAAGAATCAAGTACCTCACGTCTGTTAACTGCTTTGAAACTTTATCATTTTCTGGGA
 AAGGTGCTGGTGGTGTGCTTACGCTTTTCGGGGCCTGATGGCGTTAACTACAACCGTGCCTCCAGG
 TATGGGCCATCCCTCTGCTGCTGGTTCCTTTTCGCCTACTTAGCAGCTCACAGTTTCTATCTGTGTT
 CGAAACTGTGCTGGATACACTTTTCTGTGCTTTGCTGTGACCTGGAACAAATGATGGATCGTCGGAG
 AAACCGTACTTTATGGACCCAGGTTTTTGAGTTTTGAAAACGGACCGACCACTTAACATGCCAGGT
 CTC AAGGACACAAGAACTCACTACCAAAATGAAGAGGGGACAGAACTCCGGCCCATTTGTGAGA

ACGGGCCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGA
TTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RR214955 representing NM_001013914
 Red=Cloning site Green=Tags(s)


MGYSVVAAGRLLFGYDSFGNVCGKRNSPVEGAPLSGQDMTLKKHVFFMNACNLEVKDRGRGPMALCVS
 SCPEKQLDLEEVLQFANINGSFLCVYLSNSFNYSQSPSADRLCPRLPVPPSKPFPFLFNRCIPQTPECY
 LFASVLLNDADALHRILSGIMAGRDTILGLCAFVFALSLAMLFTFRFISTLLVHIIISLVILGLLFCVGV
 FWWLYDYDNDLSTELDTEKENMKMLAFVITVVTVLLALIFTLRKRKLVTELLRVNKAISRCP
 LLLQPLWTFAILVFFWLVWAVLLSLGTAGTAQVMEGGQVEYKPLSGIRYLWYHLIGLIWTSEFILTQC
 RMTVAGAMVTCYFNRNQNDPPARPILSSLSTLFCYHQGTAVKGSLLLTVTRIPRIVFMYIYSTVKERHSA
 WPRGEFRCSYGLWCLTNYPYHLNQDAYAAAAINGTDFCTSAKDAHTIISKNSSHLTSVNCFGNFIIFLG
 KVLVVCFSVFGGLMAFNYNRALQVWAIPLLLVAFFAYLAAHSFLSVFETVLDLFLCFVADLETNDGSS
 KPVFMDPGFLSFVKRTDHFNNARSQGHKNSLPNEEGTELRPIVR

TRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-NotI

Cloning Scheme:



ACCN: NM_001013914

ORF Size: 1812 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001013914.1](#), [NP_001013936.1](#)

RefSeq Size: 2302 bp

RefSeq ORF: 1815 bp


Locus ID: 295417

UniProt ID: [Q6AY92](#)

Cytogenetics: 2q41

MW: 67.7 kDa

Product images:



Circular map for RR214955