

Product datasheet for **RR214956L4V**

Slc16a1 (NM_012716) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Slc16a1 (NM_012716) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Slc16a1
Synonyms:	MCT1; RATMCT1; RNMCT1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_012716
ORF Size:	1482 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR214956).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_012716.2 , NP_036848.1
RefSeq Size:	2517 bp
RefSeq ORF:	1485 bp
Locus ID:	25027
UniProt ID:	P53987
Cytogenetics:	2q34
Gene Summary:	transports lactic acid and other monocarboxylic acids [RGD, Feb 2006]



[View online »](#)