

Product datasheet for RR214973L3V

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Sncg (NM_031688) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type: Lentiviral Particles

Product Name: Sncg (NM_031688) Rat Tagged ORF Clone Lentiviral Particle

Symbol: Sncg

Mammalian Cell Puromycin

Selection:

Vector:

pLenti-C-Myc-DDK-P2A-Puro (PS100092)

Tag: Myc-DDK

ACCN: NM_031688

ORF Size: 369 bp

ORF Nucleotide

The ORF insert of this clone is exactly the same as(RR214973).

Sequence:

Cytogenetics:

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

RefSeq: <u>NM 031688.1</u>, <u>NP 113876.1</u>

16p15

 RefSeq Size:
 732 bp

 RefSeq ORF:
 372 bp

 Locus ID:
 64347

 UniProt ID:
 Q63544

Gene Summary: may play a role in the neurofilament network [RGD, Feb 2006]