

Product datasheet for **RR215016L3V**

Gp1bb (NM_053930) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Gp1bb (NM_053930) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Gp1bb
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_053930
ORF Size:	618 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR215016).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_053930.4 , NP_446382.1
RefSeq Size:	947 bp
RefSeq ORF:	621 bp
Locus ID:	116727
UniProt ID:	Q9JIM7
Cytogenetics:	11q23
Gene Summary:	forms GP1b, a receptor for von Willebrand factor, when linked by a disulfide bond to platelet glycoprotein GP1alpha; may be involved in platelet aggregation [RGD, Feb 2006]



[View online »](#)