

Product datasheet for **RR215042L4V**

Scrn1 (NM_001025063) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Scrn1 (NM_001025063) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Scrn1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001025063
ORF Size:	1242 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR215042).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001025063.1 , NP_001020234.1
RefSeq Size:	1584 bp
RefSeq ORF:	1245 bp
Locus ID:	502776
UniProt ID:	Q6AY84
Cytogenetics:	4q24
Gene Summary:	Regulates exocytosis in mast cells. Increases both the extent of secretion and the sensitivity of mast cells to stimulation with calcium (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)