

## Product datasheet for **RR215074**

### Vgf (NM\_030997) Rat Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Vgf (NM_030997) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Vgf
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

**ORF Nucleotide Sequence:**

>RR215074 representing NM\_030997  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGAAACCTTACGTTGCCAGCATCCGTCCTTCTGCTTCTCTACTCATCCGGGGTTGGGAGCAG  
 CACCCCCGGGCGCTCCGATGTTTATCCTCCTCCCTCGGCTCTGAGCATAATGGGCAGGTAGCTGAGGA  
 CGCAGTGTCCCGGCCAAAGGATGACAGCGTCCCAGAGGTCGAGCGGCTCGGAATTCGAGCCTCAGGAC  
 CAGGGAGAGCTCTCCAGGGCGTGGATCCCGGGCGCTGGCCGCGTACTGTTGCAGGCACTGGACCGTC  
 CGGCTCGCCCCGGCTGTCCCGCAGGTTCCAGCAGGGAACCCGAAGAAGCAGCAGAAGCTCTGCT  
 GACCGAGTCCGTGCGCAGTCAACCCATAGCCTCCCGCATCAGAAATCCAAGCGTCCGCTGTGGCGCC  
 CCTCGCCCTCAGACTCAGGACAACGATCCCGAGGCAGACGACCGCTCAGAAGAGCTGGAGGCACTAGCAT  
 CCTTGCTCCAAGAACTTCGAGATTTTCAGTCCGAGTAATGCTAAGCGCCAGCAAGAGACGGCGGCAGCAGA  
 GACTGAAACCCGCACGCACACGCTGACCCGAGTCAATCTGGAGAGCCCCGGGCCAGAGCGGTATGGCGC  
 GCTTCCGAGGAGAGTTCCAGGCGCGCTCCCGAGCGTCTCCTCTGCCGCCCTCGGTCCCTTCTCAAT  
 TCCAGGCTCGAATGTCCGAAAACGTTCCCTTCCCGAAAACCCATCAGTTCGGGGAAGGAGTGTCTCCCC  
 TAAACACATCTTGGTGAGACTTTGACACCCCTTATCCAAGGCGTACCAAAGTCTAAGTGCCCCCTTCCCC  
 AAGGTGCGTCGGCTCGAGGGCTCATTCTGGGCGGTTCCGAGGCAGGAGAGCGCTGTTCAACAAGGT  
 TAGTCCAGGTAGAGGCAGGGAGGAGGCAGGCGGAGGCCACCCGGCAGGCCGACGCAAGAAGAGCGGCT  
 GGCCGATCTCGCTCCGACCTGCTGCTCCAGTATTTGCTGCAGGGCGGCGCCGGCAGCGCATCTGGG  
 GGTCCGGGCTGCAGGAGACGCAGCAAGAGCGGGAACGAGAGGGAGGAGGAGGCGGAGGCAGGAGCAGGAG  
 CGGTGGTGGGAGGACGAGGTGGGGAAAGAGGATGAGGAGCGGCAGAGCGGAGGCGGAGGCAGGAGGA  
 GCGGAGAGGGCGCGGCAAGACGCGCTCCTGTTCCGCGAGGAGGAGGACGGGGAAGCCGGAGCCGAGGAC  
 AAGCGCTCCAGGAGGAGGCGCCAGGCCATCGCGGAAGGATGCTGAGGGGACAGAGGAGGCGGGGAGG  
 AGGATGACGACGCAAGAGATGGATCCGACAGCATCGATAGTCTCATTGAACTGTCCACCAAACCTCCA  
 CCTGCCAGCAGCATGTGGTCAAGCATCGAAGAGGTGGAGGAGAAACGGAAGCGGAAGAAGAAGCC  
 CCTCCGAGCCGGTGCCTCCCCCAGGGCTGCCAGCCCGACCCATGTCCGCTCCCCGAGCCCCCAC  
 CTCCCGCCCGGGGATGAGTTGCCGACTGGAACGAAGTACTCCACCTGGGATCGGAGGAGGA  
 TGAGGTGTTTCCCCGGGCGCTATCACCCCTCCCAAACACATTCCGCGCGGACACTGCAGCCGCC  
 GCATCTCCCGCGCCGTCATTCCATACGCGTTGCCACCTGCGCGCCACCATCCCGATCTGGAGGCC  
 AGGCCAGGCGCGCAGGAGGAAGCGGACGCGGAGGAGCGCGGCTGCAGGAGCAGGAGGAGCTGGAGAA  
 TTACATTGAGCACGTGCTGCTGCACCGCCG

AG**CGGACCG**ACGCGTACGCGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC  
 TGGATTACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:**

>RR215074 representing NM\_030997  
 Red=Cloning site Green=Tags(s)


MKFTFLPASVLCFLLLRGLGAAPPGRSDVYPPPLGSEHNGQVAEDAVSRPKDDSVPEVRAARNSEPQD  
 QGELFQGVDPRALAALLQALDRPASPPAVPAGSQGTPEEAAEALLTESVRSQTHSLPASEIQASAVAP  
 PRPQTQNDPEADDRSELEALASLLQELRDFSPSNAKRQQTAAAEETRTHTLTRVNLSPGPERVWR  
 ASWGEFQARVPERAPLPPSVPSQFQARMSENVPLPETHQFGEVSSPKTHLGETLTPLSKAYQSLSAPFP  
 KVRRLGSLGGSEAGERLLQQGLAQVEAGRRQAEATRQAAAQEERLADLASDLLLQYLLQGGARQDLG  
 GRGLQETQQERENEREEAEQERRGGGEDEVGEEDEAAEAEAEAEERARQNALLFAEEEDGEAGAED  
 KRSQEEAPGHRKDAEGTEEGEEDDDDEMDPQTIDSLIELSTKLHLPADDVVSIIIEVEEKRRKKNNA  
 PPEVPVPPRAAPATHVRSPPPPAPARDELPDWNEVLPWDREDEVFPPGPHYFPNIRPRTLQPP  
 ASSRRRHFFHALPPARHHPDLEAQRRAQEEADAERRLQEQEELENYIEHVLLHRP

SGP**TRTRPLEQKLI**SEEDLAANDILDYKDDDDKV


**Restriction Sites:**

Sgfl-RsrII

Cloning Scheme:



Plasmid Map:



ACCN: NM\_030997

ORF Size: 1851 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_030997.1</a> , <a href="#">NP_112259.1</a>
<b>RefSeq Size:</b>	2636 bp
<b>RefSeq ORF:</b>	1854 bp
<b>Locus ID:</b>	29461
<b>UniProt ID:</b>	<a href="#">P20156</a>
<b>Cytogenetics:</b>	12q12
<b>MW:</b>	68.2 kDa
<b>Gene Summary:</b>	peptide growth factor that plays a key role in neuronal differentiation and survival; may regulate synaptic function in hippocampal neurons [RGD, Feb 2006]