

Product datasheet for **RR215085**

Sec24b (NM_001106474) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Sec24b (NM_001106474) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Sec24b
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RR215085 representing NM_001106474 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCCGATCGCC

ATGTGGGCCCGCCGGGTCCCCTACCCGGCCGGCCGGCCGGATGCCGCCAAAGCTCGGCGGAGCAG
TCTCTGGCCTGGCGCCCGCCGAGCAGAATGGTCCAGCAGAGTCAAATGCAAGTCCATCTGGGTATGG
GCTGCCCATCAAACTATATGGCCCTCAGGACTACTCTCAAGACCTGGGAAATGACCTCCTTG
CCGTTGGATAACCAAGTGTGAGAATTAATACTCTCGTCCATATACAGCACCAACACAGAATGTGGGGACTC
CCAGCTCAGCAAACCAACAGGAGCAGCTCATGTATGGCAGAGGTCCTTCTGCCCTCACATGGGGG
ATCCATGCCAGGACCTTTC AAGAGCTCCAGCATCGGGTCCATTTCGTACCGAGTGCCTCCAGCCC
TATTCCTCTCTTGGAAACCGCTACAGTAGTCTGCCACATACTCCGCCACCGTTCTGTTGCTTCTCAGG
GATACCCCTCCACTTGTAGTCATTATCCCATCTCCACTGTGTCTAATGTCGTGTATCCTAACGTTTCATA
CCCCTCGCTGCCTGCCAGTGAGCCGTATGGGCAGATGTTACCTCACAGAGTGTCTCCGCTGCCAGG
CCACTTAAAGAGTCATACTCTGGCCAAGCAGCTGTTGCCTACCCATCGAGACCTCCGCTCCACCTT
CCCAGCACCAGCAGCAGCAGCAGCAGCAGCAGCAGCAGCAGCAGCAGCAGCAGCAGCAGTATTGAGTCTCT
ACCATGGTCAGGTCAGGCTTCCACCAGCCAAAGACAGTCTCATCCGAAACAAATGGGATCCTTGGCT
ATACCGAACAGCCACCCAGCAATTAACGAAATGTCCAACCACCTAAGCCAGCTCCGTAAGTACGACAG
TTTTGCCAGGACCTCATCAACAAGATGCCACCTGCCCAAGTACCCAGTTGGGCTGTGCCGTCTGC
TCCACCCCCCGCAGGAGATGCAGACAAAAGGCATGCAGTATGGTGACTATGTTAATAATCAAGTACG
TCCACACCAACTCCCTTGTATCAGCTTCCGATGATGAGGAAGAGGAGGATGAGGATGAGGAAGCAGGTG
TTGACAGCTCTTCAACAAGCAGCGCTTCGCCCTTGCCCAACAGTTATGATGCCTTGGAAAGGAGGA
CTACCCAGATATGCATTCTTTCAGCAAGCAGCCCGTTCTGATCGTGCCCGGAAACCTAACTCCACC
CTCGTCCCCTCCACTGCTGCTCAGCAGCCAGGTAGCAAAGCTTTGGCTATGGTACCCAACG
TGCAGCCTGTTATCAGAATGCAGCAGCACCACCACAGCCACCCAGCCAGGGGCGCAGGTATTCTGG
ATATCCTCAGCAGTATCCAGGTGTGCACCAACTGCTCAGGTCTGGGAGGGCTAAGCCTGCAGAGCTCT
CCACAGCCAGAAAAGCCTACGGCCTGTGAACCTCACCCAGGAGAAGAACATCTTACCTGCAACTCCCATT
GGGCACCTGTACCTAACCTGAGTGCAGAGCTCAGCAAAGTTAAACTGTAGCCCAGACTCGTTTCGATGTAC
CTTGACCAGTATCCCACAGACAGGCTTTACTGAATAAAGCTAAGCTTCCTTGGGACTGCTGCTGCAC



[View online »](#)

CCCTTCAGAGACCTGACGCAATTGCCAGTGATAACCTCCAACACCATTGTGAGGTGCCGGTCTGTGCGGACGTACATCAACCCTTTTGTGTCCTTCATCGATCAGCGCAGGTGGAAATGCAACCTGTGCTACAGAGTGAA
TGATGTTCCAGAAGAGTTTTTGTACAACCCCTGACTCGGTCTTACGGAGAGCCTCATAAACGGCCGGAA
GTTCAGAATCAACTGTGGAGTTCATCGCTTCCCTCAGACTACATGCTCCGGCCTCCTCAGCCTGCAGTCT
ACTTGTGTTTTAGATGTGTCTATAACGCAGTGGAAAGCTGGTACTTGACAGTTCTATGCCAGTCCTT
ATTAGAAAACCTAGACAAGCTTCCCTGGAGACTCACGAACAAGAATAGGATTATGACTTTTCGACAGCACT
ATTCATTTCTACAATTTACAAGAAGGATTATCCCAGCCTCAGATGTTGATTGTATCTGATATAGACGATG
TTTTCTTGCCACACCGGATAGTCTGCTTGTGAATCTGTATGAAAGTAAAGAGCTTATAAAAGACTTACT
GAACGCATTACCAAATATGTTTCATCAATACCAGAGAAAACGCACAGTGCTCTTGGACCTGCACTTCAGGCT
GCCTTCAAACCTAATGTCTCCACGGGCGGACGAGTGTCTGTGTTTCAGACACAGTTACCTTCCTTGGGGG
CAGGACTTCTGCAGTCCAGAGAAGATCCGAATCAGAGGTCGAGCACCAAGGTGGTGCACCATCTCGGCC
TGCCACCGACTTTTATAAGAACTTGCTTTAGATTGCTCAGGACAGCAGACCCGAGTGGATTTGTCCTG
CTAAGTTCACAGTATTCTGACCTCGTTCTCTGGCCTGCATGTCCAAGTACTCCGACGGTGCATATTCT
ACTACCATCGTCCACTCCACCACAACCTTCCGAAGCAGAGAAGTTACAGAAAGACCTAAAACGGTA
CCTCACCAGGAAGATCGGATTTGAAGCTGTTATGAGGATAAGATGTACCAAAGGACTTTTCGATGCACACT
TTCCACGGGAACCTCTTTGTCCGGTCCACTGACTTGCTGTCTCTTGCCAACATCAATCCTGATGCTGGGT
TTGCAGTGCAGTTGTCAATTGAGGAAAGCCTGGCGGACACTTCCCTGGTGTGCTTCCAGACTGCCCTGCT
GTACACATCAAGCAAAGGTGAGCGGAGAATTCAGTTCACACACTGTGCTTGCCAGTGGTGGAGTTCACTG
GCTGATGTTTATGCAGGAGTGGATGTGCAAGCTGCCATCTGCCTTCTGGCCAACATGGCTGTGGACCGCT
CAGTGTATCCAGCCTGTCCGATGCCAGGGACGCCCTGGTGAATGCCGTGGTGGATCCCTGTCTGCCTA
CAGCTCTGCAGTGTGAGTGTGCCCCGCTCCACCCTGACAGCACCGAGCTCACTGAAGTGTCTCCCTCTC
TATGTTCTGGCACTTCTCAAACAGAAAGCATTAGAACAGGTACAAGCACGCGCCTGGACGACCGTGTGT
ACGCCATGTGCCAGATGAAGTCGCAGCCACTTGTGCATCTAATGAAAATGATTTCACCAACTTGTACAG
GATAGACCGGCTGACGGATGAGGGTGCCATACATGTCAACGACAGGGTAGTACCCAGCCACCTCTTCAA
AAACTGTCTGCAGAGAACTGACCAGAGAAGGTGCTTTCCTCATGGACTGTGGCTCAGTTTTCTACATTT
GGATTGGCAAAGGCTGTGACAGTAATTTTCATAGAGAATGTGCTTGGGTATCCTGATTTTGGATCAATACC
ACAGAAAATGACACATCTTCCAGAGCTAGACACACTTCCGTCAGAGAGAACTAGATCCTTTGTAACCTGG
CTCAGAGACAGCCGACCTTAAGCCCTGTCTTACGTGGTAAAAGATGAGAGTCTGCCAGAACAGACT
TTTTCCAGCATTGGTGGAAAGACCCGACCGAAGCCGCACTCTCCTACTATGAGTTTCTGATTACGTTCA
GCAGCAGGTTTGAAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA


Protein Sequence: >RR215085 representing NM_001106474
 Red=Cloning site Green=Tags(s)


MSAPAGSPHPAAGARMPPKLGGA V SGLAPPQQNGPAQSQM V P S G Y G L P H Q N Y M A P S G H Y S Q G P G K M T S L
 P L D N Q C E N Y Y S R P Y T A P T Q N V G T P S S A N Q P G A Q L M Y G R G P S A P H M G A S M P G P F Q G A P A S A S H S Y P S A S Q P
 Y S S L G N R Y S S P A T Y S A T A S V A S Q G Y P S T C S H Y P I S T V S N V V Y P N V S Y S L P A S E P Y Q M F T S Q S A P P P A R
 P L K E S Y S G P S T A V A Y P S R P P P P S Q H Q Q Q Q Q Q Q Q Q S H S G Y S S L P W S G P G L P P A Q D S L I R N Q M G S L A
 I P N S H P A I N E N V Q P P K P S S V V A T V L P G P S S T R M P P A P S H P V G P V S A P P P P E Q M Q T K G M Q Y G D Y V N N Q A S
 S T P T P L S S A S D D E E E E D E E A G V D S S T T S S A S P L P N S Y D A L E G G S Y P D M H S S A S S P V P D R A P E P N S T
 L V P T P T A A Q P A K V A K P F G Y G Y P T L Q P A Y Q N A A P P T T A H P S G P A Y S G Y P Q Q Y P G V H Q L S S G L G L S L Q S S
 P Q P E S L R P V N L T Q E K N I L P A T P I W A P V P N L S A E L S K L N C S P D S F R C T L T S I P Q T Q A L L N K A K L P L G L L L H
 P F R D L T Q L P V I T S N T I V R C R S C R T Y I N P F V S F I D Q R R W K C N L C Y R V N D V P E E F L Y N P L T R S Y G E P H K R P E
 V Q N S T V E F I A S S D Y M L R P P Q P A V Y L F V L D V S H N A V E A G Y L T V L C Q S L L E N L D K L P G D S R T R I G F M T F D S T
 I H F Y N L Q E G L S Q P Q M L I V S D I D D V F L P T P D S L L V N L Y E S K E L I K D L L N A L P N M F I N T R E T H S A L G P A L Q A
 A F K L M S P T G G R V S V F Q T Q L P S L G A G L L Q S R E D P N Q R S S T K V V H H L G P A T D F Y K K L A L D C S G Q Q T A V D L F L
 L S S Q Y S D L A S L A C M S K Y S A G C I F Y Y P S F H S T H N P S Q A E K L Q K D L K R Y L T R K I G F E A V M R I R C T K G L S M H T
 F H G N F F V R S T D L L S L A N I N P D A G F A V Q L S I E E S L A D T S L V C F Q T A L L Y T S S K G E R R I R V H T L C L P V V S S L
 A D V Y A G V D V Q A A I C L L A N M A V D R S V S S S L S D A R D A L V N A V D P L S A Y S S A V S S V P R S T L T A P S S L K L L P L
 Y V L A L L K Q K A F R T G T S T R L D D R V Y A M C Q M K S Q P L V H L M K M I H P N L Y R I D R L T D E G A I H V N D R V V P Q P P L Q
 K L S A E K L T R E G A F L M D C G S V F Y I W I G K G C D S N F I E N V L G Y P D F G S I P Q K M T H L P E L D T L P S E R T R S F V T W
 L R D S R L S P V L H V V K D E S P A R T D F F Q H L V E D R T E A A L S Y Y E F L I H V Q Q V C K

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:


ACCN: NM_001106474

ORF Size: 3726 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001106474.1](#), [NP_001099944.1](#)

RefSeq Size: 5001 bp

RefSeq ORF: 3729 bp

Locus ID: 295461

Cytogenetics: 2q42

MW: 134.7 kDa