

Product datasheet for **RR215087**

Rbpj (NM_001106631) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Rbpj (NM_001106631) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Rbpj
Synonyms:	Rbpj
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>RR215087 representing NM_001106631
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGCGAAATTATTTAAAAGAACGAGGGGATCAAACAGTACTATTCTTCATGCAAAAAGTGCACAGAAGT
 CATATGGAAATGAAAAACGATTTTTTGCCTCCTCCTTGTGTATCTTATGGGACAGTGGTTGGAAGAA
 AAAAAAGAACAATGGAACGAGATGGTTGTTCTGAACAAGAGTCTCAACCCTGTGCATTTATTGGAATA
 GGAATAGTGACCAAGAAATGCAGCAGCTGAACCTGGAAGGGAAGAACTACTGCACAGCCAAAACATTGT
 ACATATCTGATTAGACAAGAGGAAACATTTTCATGTTGTCTGTAAGATGTTCTATGGCAACAGTATGA
 CATCGGTGTTCCTTAGCAAGCGGATAAAGGTCATCTCAAACCTTCCAAAAAGAAGCAATCATTGAAA
 AATGCTGACTTGTGCATTGCCTCAGGGACAAAGGTGGCGTGTAAATCGACTTCGGTCCCAGACCGTTA
 GCACCAGATACCTGCACGTAGAAGGAGGAATTTCCATGCTAGTTCACAGCAGTGGGAGCATTTCAT
 TCATCTCTGGATGATGATGAATCGGAAGGAGAGGAATTCACAGTTAGAGATGGTACATCCATTATGGA
 CAGACTGTCAAGCTGGTGTGCTCAGTACTGGCATGGCACTCCCAAGACTGATAATTAGGAAAGTTGATA
 AGCAGACAGCATTACTGGATGCAGACGACCCTGTATCACAACTCCACAAATGTGCATTTTACCTTAAGGA
 CACAGAAAGAAATGACTTGTGCCTTTCTCAAGAAAGAATAATCCAATTTCAAGCCACTCCATGTCCAAAA
 GAACAAAATAAGAAATGATAAATGATGGAGCGTCCTGGACAATCATCAGCACAGACAAGGCTGAGTACA
 CGTTCTATGAGGGAATGGGCCCTGTCTTGCCCGGTCACTCCTGTGCCTGTCTAGAAAGTCTCCAGTT
 GAATGGCGCGGGGATGTAGCCATGCTTGAACCTCACAGGACAAACTTTACTCCAATTTACGAGTGTGG
 TTTGGGGATGTAGAGGCCGAAACAATGTACAGATGTGGAGAGAGCATGCTCTGTGTGGTCCCAGACATCT
 CTGCATCCGAGAAGGTTGGCGATGGTCCGGCAGCCAGTCCAGGTTCCAGTAACTTGGTCCGTAATGA
 TGGAAATCATTACTCCACAGCCTTACCTTTACCTACACCCAGAACAGGGCCAAGGCCGCACTGCAGC
 GCAGCGGAGCCATTCTCAGAGCCAACCAAGCCAAGTGCCACCAATGAGTCAAACACAAAACAGCGAGG
 GAAATTACACAAATGCCAGCACAAATCTACCAGTGTCACGTCGTCCACAGCAACCGTGGTGTCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RR215087 representing NM_001106631
 Red=Cloning site Green=Tags(s)


MRNYLKERGDQTVLILHAKVAQKSYGNEKRFFCPPPCVYLMGSGWKKKKEQMERDGCSEQESQPCAFIGI
 GNSDQEMQQLNLEGNKYCTAKTLYISDSKDRKHFMLSVKMFYGNSSDDIGVFLSKRIKVISKPSKKKQSLK
 NADLCIASGKVALFNRLRSQTVSTRYLHVEGGNFHASSQWGFYIHLDDDESEGEFTVRDGYIHYG
 QTVKLVCSVTGMALPRLIIRKVDKQTALLDADDPVSQLHKCAFYLKDERMYLCLSQERIIQFQATPCPK
 EQNKEMINDGASWTIISTDKAEYTFYEGMGPVLPVTPVPVVEVSLQLNGGGDVAMLELTGQNFNPNLRVW
 FGDVEAETMYRCGESMLCVVPDISAFREGWRWVRQPVQVPVTLVRNDGIIYSTSLTFTYTPPEGPRPHCS
 AAGAILRANSSQVPTNESNTNSEGNYTNASTNSTSVTSSTATVVS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:

Plasmid Map:


ACCN: NM_001106631

ORF Size: 1395 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001106631.1](#), [NP_001100101.1](#)

RefSeq Size: 4279 bp

RefSeq ORF: 1398 bp

Locus ID: 679028

Cytogenetics: 14q11

MW: 52 kDa