

## Product datasheet for **RR215091**

### **Clptm1 (NM\_001106232) Rat Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	Clptm1 (NM_001106232) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Clptm1
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

**ORF Nucleotide Sequence:**

>RR215091 representing NM\_001106232  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**GCGATCGCC**

ATGGCGGGCGCAGGAGCGGACGGGGCCGGCAGCGCCGTGGTGGCGGGGGGAGGCAGCTCTGGTC  
 AGGTGACCAGCAATGGCAGCATTGGGAGGGACACACCAGCCGAGACCCAGCCTCAGAACCCCTCCACCCCA  
 GCCAGCCCCAACGCCTGGCAGGTCATCAAAGGGTTCTGTTTCAGGATCTTCATCATCTGGGCCATCAGC  
 AGCTGGTTCGCCGAGGGCCATCCCCTCAGGACCAGTCAGGCCCTGGAGGAGCCCCACGTGTCGCCAGCC  
 GAACTTATTCCCAAAGACTCTAATGAACCTGCACGTGTACATCTCAGAGCAGGAGCACTTCACAGA  
 CTTCAATGCCACCTCGGCGCTCTTCTGGGAGCAACATGACCTGGTCTATGGTACTGGACCAGTGGAGAG  
 AACTCCGATGGTGTACGAGCACTTTCAGAGCTTGACATCCCACAGAGTGTACAGCAGAATGGCTCCA  
 TCTATATCCATGTCTATTTACCAAGAGTGGCTTCCACCCAGACCCCGCAGAAGGCCCTGTACCGCCG  
 GCTTGCCACCGTCCACATGTCCCGGATGATCAACAAATACAAGCGTAGACGATTTAGAAAACCAAGAAC  
 CTAAGTACAGGAGAGACAGAAGCTGACCCGAAATGATCAAGAGGGCTGAGGACTATGGCCCTGTAGAAG  
 TGATCTCACACTGGCACCCCAACATCACCATCAATATCGTGGATGACCACACACCATGGGTGAAGGGCAG  
 CGTGCCACCCCACTGGACCAATATGTGAAGTTTGTGCGGTGAGTGGTACTACTACCCATCATCTAT  
 TTCAATGACTACTGGAACCTGCAGAAGGATTACTACCCATTAATGAGAGCCTAGCCAGTCTACCACTGC  
 GCGTCTCTTTCTGCCCTGTCACTGTGGCGCTGGCACTTTATGCTGCACAGAGTACCAAGTCAACCTG  
 GAATTTCTTAGGTGATGAGCTGTATGAGCAATCGGACGAGGAACAGGATTCAGTGAAGTTGCCCTGCTG  
 GAGACCAGCCCTATCTTCTGGCACTTACCATCATCGTCTCCATTGTCCACAGTGTCTTCGAGTTCTGCT  
 CTTTCAAAAATGACATTCAGTTCTGGAACAGTCGTCAGTCTCTGGAGGGGCTGTCTGTGCGGTCTGTCTT  
 CTTTGGCGTTTTCCAGTCCTTTGTGGTTCTCCTTTACATCCTGGACAATGAGACCAACTTCGTGGTACAG  
 GTCAGCGTCTTCATTGGGGTCTCATTGACCTCTGGAAGATTACCAAGGTCATGGATGTCGGCTGGACC  
 GGGAGCACAGAGTGGCAGGGATCTTCCCTGTCCGACATTAGAGACAAGTCCACCTACATTGAGTCTTC  
 AACCAAGTATATGATGACATGGCGTTCGGTACCTGTGATGGATCCTCTTTCCCTGCTAGGCTGCTAT  
 GCCGTCTACAGCCTTCTGTACCTGGAGCACAAGGGCTGGTACTCCTGGGTGCTCAGCATGCTCTACGGCT  
 TCCTGCTGACCTTTGGCTTCATCACCATGACACCCCACTGTTTCACTAATAAGCTCAAGTCTGTGGC  
 CCACTTGGCATGGCGAATGCTCACCTACAAGGCCCTCAACTTTTATTGATGACTTGTTCCTTTGCTTGTG  
 ATCAAGATGCCTGTGATGATACCGGATTGGCTGCCTGCGGGACGATGTGGTTTTTCTCATCTACCTTACC  
 AACGGTGGATCTACCGCGTCGACCCACCCGTGTGAACGAGTTCCGGTATGAGTGGAGAGGATGTCTCAGC  
 AGCGGCATCCGGGGCCAGGCATCCACTGCAGCAGGGGCTCTCACACCTTACCCAGCACAGCTGTACCC  
 GGAGAGGATGCCTCCACAGTCCCAACGGCTACTCCTGGGGCTGCACGGCCAGCCAACCCAGGAAGACC  
 CTCAAAACCAGCAGAGGACAAGAAAAGGAT

**ACGCGT**ACGCGGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR215091 representing NM\_001106232  
 Red=Cloning site Green=Tags(s)


MAAAQEADGAGSAVVAAGGGSSGQVTSNGSIGRDTPAETQPQNPPQPAPNAWQVIKGVLFRIFIWAIS  
 SWFRRGSPQDQSGPGGAPRVASRNLFPKDTLMNLHVYI SEHEHFTDFNATSALFWEQHDLVYGDWTSGE  
 NSDGCYEHFAELDIPQSVQONGSIYIHVYFTKSGFHPDPRQKAL YRRLATVHMSRMINKYKRRRFQKTKN  
 LLTGETEADPEMIKRAEDYGPVEVISHWHPNITINIVDDHTPWVKGSPVPPPLDQYVKFDAVSGDYPIIY  
 FNDYWNLQKDYYPINESLASLPLRVSFCLSLWRWQLYAAQSTKSPWNFLGDELYEQSDEEQDSVKVALL  
 ETPYLLALTIIIVSIVHSVFEFLAFKNDIQFWSRQSLLEGLSVRSVFFGVFQSFVLLYILDNETNFVVQ  
 VSVFIGVLIDLWKITKMDVRLDREHRVAGIFPCPTFRDKSTYIESSTKVYDDMAFRYLSWILFPLLGCY  
 AVYLLYLEHKGWYSWLSMLYGFLLTFGFITMTPQLFINYKLSVAHLPWRMLTYKALNTFIDDLFAFV  
 IKMPVMYRIGCLRDDVFFIYLQYRWIYRVDPTRVNEFGMSGEDVSAASGAQASTAAGALTSPSPSTAVT  
 GEDASTVPTATPGACTASQPQEDPPKPAEDKKK

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV


Restriction Sites:  
 Cloning Scheme:

SgfI-MluI

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF

**Plasmid Map:**


**ACCN:** NM\_001106232

**ORF Size:** 1992 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_001106232.1](#), [NP\\_001099702.1](#)  
RefSeq Size: 2453 bp  
RefSeq ORF: 1995 bp  
Locus ID: 292696  
Cytogenetics: 1q21  
MW: 75.3 kDa