

Product datasheet for **RR215099**

Serpina12 (NM_138825) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Serpina12 (NM_138825) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Serpina12
Synonyms: OI-64; OI64
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR215099 representing NM_138825
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAACCTCGTGCTGGGCTGGGTCTGTTTCTGGCAGGCCTCCTCACTGTGAAAGGTCTCCTGCAGGACA
GAGATGCTCCAGATACATATGAGTCTCCAGTCAGGTCCTCAAGAGTGGAGGGCAAGAAGGATGCCGGGA
GCTTACTCGGCACAACATGGAATTTGGCTTCAAGCTACTGCAGAGGCTGGCCAGCAACAGCCGTCAGGGG
AACATCTTCTTATCTCCACTGAGCATCTCAACTGCCTTCTCCATGCTGTCTCTGGGGGCCAAAACTCCA
CACTGGAGGAAATCAGGGAGGGTTTCAATTTCAAGGAGATGTCAGACAGGGACATGCACATGGGTTTCCA
CTACCTTCTCCAAAAGCTCAACCGGGAGACACAGGACGTAAGAGATGAGCATAGGCAATGCTTTATTTATG
GACCAGAGGCTGAGGCCCAACAGAGGTTTCTGAAACTGGCTAAGAACTTGATGATGCGGACATGATAC
TCACTAATTTCCAGGACTTAGAAAACACTCAGAAGAATATCAACAAGTACATCAGTCGGAAAACCCACAA
CAGAATCGAAAACATGGTCAAAAACATAGACCCAGGGACGGTGATGCTCCTCACAACTATATTTACTTT
CAAGGCAGGTGGCAGTATGAATTCGACCCAAAACAAACAAAGAGGAGGACTTCTTCATAGAGGAAGGCA
AGACCGTGAAGGTGCCATGATGTTTCAAAGGGGCATGTACGACATGGCCTATGACAGTCAGCTCTCCTG
CACCATCTTGAAATGCCCTACCGAAGAACATCACGGCCACCTTTGTCCTTCTGACAGTGCCAACTG
AGGCTCCTGGAGCAGGGCTTGAAGCAGACATCTTTGCCAAATGAAAAGCCTTCTCTCTAAAAGAGTGG
TAGATGTGTGGGTGCCAGGCTCCACATCTCGGCTACTTACAACATGAAGAAGGTTCTGTCCCGCTGGG
CATCAGCAAAATCTTTGAAGAGCATGGTGACCTCACCCGCATCTCGTCTCATCGAAGCCTGAAGGTGGGA
GAGGCTGTGCACAAGGCTGAGCTGAGGATGAATGAGAAGGGCACAGAAGGAGCTGCAGGCTCGGGAGCCC
AGACACTGCCATGGAGACCCACGAAGGATGAAGCTAAACGCACCCTTCTGATGATGATTTATGAGAA
CCTCATGCCTTCCATGATCTTCTTGCAAGGATTTATAACCCTAGTGAGAAA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR215099 representing NM_138825
 Red=Cloning site Green=Tags(s)


MNLVLGLGLFLAGLLTVKGLLQDRDAPDITYE SPVRVQEWGKKDARELTRHNMEFGFKLLQRLASNSRQG
 NIFLSPLSISTAFSMLS LGAQNSTLEEIREGFNFKEMSDRDMHMGFHYLLQKLNRETQDVKMSIGNALFM
 DQRLRPQQRFLKLAKNL YDADMILTNFQDLENTQKNINKYISRKTHNRIENMVKNIDPGTVMLLTNYIYF
 QGRWQYEFDPKQTKEDFFIEEGKTVKVPMMFQRGMYDMAYDSQLSCTILEMPYRRNITATFVLPDSGKL
 RLLEQGLQADIFAKWKSLLSKRVVDVWVPRLHISATYNMKKVLSRLGISKIFEEHGLTRISSHRSLKVG
 EAVHKAELRMNEKGTEGAAGSGAQLPMETPRRMKNAPFLMMIYENLMPSMIFLARIYNPSEK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_138825

ORF Size:	1242 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_138825.2 , NP_620180.2
RefSeq Size:	1487 bp
RefSeq ORF:	1245 bp
Locus ID:	191570
Cytogenetics:	6q32
MW:	48 kDa
Gene Summary:	human homolog is a serine protease inhibitor [RGD, Feb 2006]