

Product datasheet for **RR215100L3V**

PPP3R2 (NM_021701) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	PPP3R2 (NM_021701) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	PPP3R2
Synonyms:	Cblp
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_021701
ORF Size:	528 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR215100).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_021701.2 , NP_067733.1
RefSeq Size:	739 bp
RefSeq ORF:	531 bp
Locus ID:	29749
UniProt ID:	P28470
Cytogenetics:	5q22
Gene Summary:	may act as a testis-specific calcium dependent protein serine/threonine phosphatase isozyme [RGD, Feb 2006]



[View online »](#)