

Product datasheet for **RR215112**

Gpr65 (NM_001106751) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Gpr65 (NM_001106751) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Gpr65
Synonyms: Tdag8
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR215112 representing NM_001106751
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGACGATGAACAGCACATGTGTTGAAGAACAGCATGACCTGGACCACTATTTGTTCCAGTGGTCTACA
TATTTGTGTTTATAGTCAGCGTCCAGCCAACATCGGCTCTTTATGTGTATCTTTCTGCAAGCGAAGAA
GGAAAATGAACTAGGAGTTTACCTCTTCAGTCTGTCACTGTCAGATCTGCTGTATGCGTTGACGCTGCCT
CTCTGGATAAACTACACTTGGAAATAAGACAACCTGGACGCTCTCTCCACCTTGTGCAAAGGAAGGTTT
TCTTCACCTACCTGAACTTTTACAGCAGCACGGCCTTCCTCACTTGCATTGCCTTGGACCGCTATTTAGC
AGTCGTCTACCCTCTAAGGTTTTCTTTCTAAGAACAAGAAGATTTGCATTTATCACCAGCCTGTCTATC
TGGATATTGGAATCCGCTTTAACTCCATCCTTCTGTGGGAAGACGAAACAAGTGTGCAATACTGTGATG
CAGAGAAATCTAATTTCACTCTCTGCTATGACAAATACCCTCTGGAGAAGTGGCAGATAAACCTTAACCT
GTTCCAGGACTTTCCTGGGCTACGCAATACCCTTGGTCACCATCATGATCTGCAACCATAAAGTCTACCGC
GCCGTGAAGCACAACCAAGCTACGAAAAACAAGAAAGGAGACGGATCATAAAGCTGCTTCCAGCATCA
CGCTGACTTTCATCTTATGCTTTACCCCTTCCACGTGATGGTGTCTATCCGCTGCATTTTAGAGCGCAA
CGTGAACCTCAATAACAAGACTGGAAAGCAGACATTTACGGTGTACAGAATCACAGTGGCCCTGACGAGT
CTAAACTGTGTCGCCGATCCAATTCTCTACTGTTTCGTGACGAAACCGGGAGATCCGATATGTGGAACA
TATTAAGTGTGTGCTAGGAAGCGCAATAGATCCCGAGGACAAAAAAGAGACATACTTTCTGTGCCAC
AAGAGATACCATAGAAGTGGAGATTATAGAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR215112 representing NM_001106751
 Red=Cloning site Green=Tags(s)

MTMNSTCVEEQHDLDHLYLFPVVYIFVFIIVSVPANIGSLCVSFLQAKKENELGVYLFSLSLDLLYALTLP
 LWINYTWNKDNWTLSPTLCKGSVFFTYLNFYSSAFLTCIALDRYLAVVYPLRFSFLRTRRFAFITSLSI
 WILES VFNSILLWEDETSVEYCDAEKSNFTLCYDKYPLEKWQINLNLFRFLGYAIPLVTIMICNHKVYR
 AVKHNQATENNERRRIKLLASITLTFILCFTPFHVMVLIRCIERNVNFNNKTGKQTFVTYRITVALTS
 LNCVADPILYCFVTETGRSDMWNILKVCARKNRNRSRGQKRDILSVSTRDITIELEIIE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001106751

ORF Size: 1011 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001106751.2 , NP_001100221.1
RefSeq Size:	1014 bp
RefSeq ORF:	1014 bp
Locus ID:	299242
Cytogenetics:	6q32
MW:	39.3 kDa