

Product datasheet for **RR215116L3V**

Itpa (NM_001107774) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Itpa (NM_001107774) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Itpa
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001107774
ORF Size:	594 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR215116).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001107774.1 , NP_001101244.1
RefSeq Size:	1184 bp
RefSeq ORF:	597 bp
Locus ID:	311422
UniProt ID:	D3ZW55
Cytogenetics:	3q36
Gene Summary:	human homolog is a nucleoside triphosphate pyrophosphatase; hydrolyzes inosine triphosphate and deoxyinosine triphosphate to the monophosphate nucleotide and diphosphate [RGD, Feb 2006]



[View online »](#)