

## Product datasheet for RR215125L4V

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## **Gpr26 (NM 138841) Rat Tagged ORF Clone Lentiviral Particle**

**Product data:** 

**Product Type: Lentiviral Particles** 

**Product Name:** Gpr26 (NM 138841) Rat Tagged ORF Clone Lentiviral Particle

Symbol: Gpr26

**Mammalian Cell** Puromycin

Selection:

Vector:

pLenti-C-mGFP-P2A-Puro (PS100093)

mGFP Tag:

ACCN: NM\_138841 **ORF Size:** 1011 bp

**ORF Nucleotide** 

OTI Disclaimer:

The ORF insert of this clone is exactly the same as(RR215125).

Sequence:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

RefSeq: NM 138841.1, NP 620196.1

1q41

RefSeq Size: 1172 bp RefSeq ORF: 1014 bp Locus ID: 192153 **UniProt ID:** <u>Q9Q</u>XI3 Cytogenetics:

**Gene Summary:** orphan G protein-coupled receptor that may be involved in regulating cell growth and shape

[RGD, Feb 2006]

