

Product datasheet for RR215136L4

Lipg (NM_001012741) Rat Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Lipg (NM_001012741) Rat Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Lipg
Synonyms:	lipase
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR215136).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.



[View online »](#)

Locus ID: 291437

UniProt ID: [Q8VBX1](#)

Cytogenetics: 18q12.2

Gene Summary: Has phospholipase and triglyceride lipase activities. Hydrolyzes high density lipoproteins (HDL) more efficiently than other lipoproteins. Binds heparin (By similarity).[UniProtKB/Swiss-Prot Function]