

## Product datasheet for **RR215187L3V**

### Ube2i (NM\_013050) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Ube2i (NM_013050) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Ube2i
Synonyms:	UbcE2A
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_013050
ORF Size:	474 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR215187).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_013050.1</a> , <a href="#">NP_037182.1</a>
RefSeq Size:	1092 bp
RefSeq ORF:	477 bp
Locus ID:	25573
UniProt ID:	<a href="#">P63281</a>
Cytogenetics:	10q12
Gene Summary:	ubiquitin-conjugating enzyme (E2 class), catalyzes the attachment of ubiquitin-like protein Sumo1 to other proteins; may be involved in the regulation of cell growth and differentiation [RGD, Feb 2006]



[View online »](#)