

## Product datasheet for **RR215214L4V**

### Sarnp (NM\_001033070) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

|                           |  |
|---------------------------|--|
| Product Type:             | Lentiviral Particles   |
| Product Name:             | Sarnp (NM_001033070) Rat Tagged ORF Clone Lentiviral Particle  |
| Symbol:                   | Sarnp  |
| Synonyms:                 | Cip29; RGD1305692  |
| Mammalian Cell Selection: | Puromycin  |
| Vector:                   | pLenti-C-mGFP-P2A-Puro (PS100093)  |
| Tag:                      | mGFP   |
| ACCN:                     | NM_001033070   |
| ORF Size:                 | 630 bp   |
| ORF Nucleotide Sequence:  | The ORF insert of this clone is exactly the same as(RR215214).   |
| OTI Disclaimer:           | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a> |
| OTI Annotation:           | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.   |
| RefSeq:                   | <a href="#">NM_001033070.1</a> , <a href="#">NP_001028242.1</a>  |
| RefSeq Size:              | 1244 bp  |
| RefSeq ORF:               | 633 bp   |
| Locus ID:                 | 362819   |
| UniProt ID:               | <a href="#">Q498U4</a>   |
| Cytogenetics:             | 7q11   |



[View online »](#)

**Gene Summary:**

Binds both single-stranded and double-stranded DNA with higher affinity for the single-stranded form. Specifically binds to scaffold/matrix attachment region DNA. Also binds single-stranded RNA. Enhances RNA unwinding activity of DDX39A. May participate in important transcriptional or translational control of cell growth, metabolism and carcinogenesis. Component of the TREX complex which is thought to couple mRNA transcription, processing and nuclear export, and specifically associates with spliced mRNA and not with unspliced pre-mRNA. TREX is recruited to spliced mRNAs by a transcription-independent mechanism, binds to mRNA upstream of the exon-junction complex (EJC) and is recruited in a splicing- and cap-dependent manner to a region near the 5' end of the mRNA where it functions in mRNA export to the cytoplasm via the TAP/NFX1 pathway (By similarity).[UniProtKB/Swiss-Prot Function]