

Product datasheet for **RR215265L4V**

Pigr (NM_012723) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pigr (NM_012723) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Pigr
Synonyms:	pIgA-R; RNPIGR2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_012723
ORF Size:	2307 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR215265).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_012723.4 , NP_036855.3
RefSeq Size:	3929 bp
RefSeq ORF:	2310 bp
Locus ID:	25046
UniProt ID:	P15083
Cytogenetics:	13q13
Gene Summary:	contains a trinucleotide repeat-rich fragment in the 3' untranslated region that may regulate gene expression [RGD, Feb 2006]



[View online »](#)