

Product datasheet for **RR215274**

Fstl3 (NM_053629) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Fstl3 (NM_053629) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Fstl3
Synonyms: Flrg
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR215274 representing NM_053629
Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGCGTCTGGGGCACTGTGGCCGCTGCTTTGGGGAGCCCTGGTCTGGGCGGTGGGATCCGTGGGCGCCG
TGATGGGCTCGGGGATTCTGTGCCGGTGGTGTGTGCTGGCTCCAGCAGGGCAAGGAGGCCACCTGCAG
TCTGGTGTGAAGACACAGGTCAAGCCGGGAGGAGTGTGCGCTTCCGGCAACATCAACACAGCCTGGTCC
AACTTCACCCACCCGGGCAATAAAATCAGCCTGCTAGGATTCCTGGGCCTTGCCACTGCCTACCATGCA
AAGATTCTTGCAGCGAGTGGAGTGCAGCCCGGCAAGGCGTGCAGCATGTTGGGGGGCCGTCCACACTG
CGAGTGCCTGTCCAACCTGCGAGGGGGTCCCGCGGGCTCCAGGTCTGCGGCTCTGATGGCGCCACCTAC
CGGGACGAATGCGAACTGCGCACCCGCGCTGTGCGGACACCCAGACCTGCGCGTCATGTACCGCGGCC
GCTGCCAAAAGTCTGCGCTCAGGTAGTGTGCCCGCGCCCGAGTCGTGCCTTGTGGATCAGACTGGCAG
CGCACACTGCGTGGTGTGTCGCGCTGCACCCTGCCAGTACCTCCCAACCCTGGCCAAGAAGTCTGCGGC
ACAACAACGTTACCTACATCTCGTGTGTCACCTGCGCCAGGCCACTTGCTTTCTGGGCGTCCATTG
GAGTTCGGCACCCGGGCATCTGCACAGGTGGCCCAAAGTACCAGCAGAGGAGGAAGAGAAGTCTCGTG

ACGCGTACGCGGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)


Protein Sequence: >RR215274 representing NM_053629
Red=Cloning site Green=Tags(s)

MRPGALWPLLWGALVWAVGSVGAVMGS GDSPGGVCWLQQGKEATCSLVLKTQVSRECCASGNINTAWS
 NFTHPGNKISLLGFLGLVHCLPCKDSCDGECPGPKACRMLGGRPHCECVSNCEGVPAGFQVCGSDGATY
 RDECELRTARCRGHPDLRVMYRGRCKSKAQVVCPRPQSLVDQTGSAHCVV CRAAPCPVPPNPGQELCG
 NNNVTYISSCHLRQATCFLGRSIGVRHPGICTGGPKVPAEEEEENFV


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_053629
ORF Size: 768 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_053629.3 , NP_446081.1
RefSeq Size:	771 bp
RefSeq ORF:	771 bp
Locus ID:	114031
UniProt ID:	Q99PW7
Cytogenetics:	7q11
MW:	27.1 kDa
Gene Summary:	may have a role in spermatogenesis [RGD, Feb 2006]