

Product datasheet for **RR215276**

Kcnab2 (NM_017304) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Kcnab2 (NM_017304) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Kcnab2
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR215276 representing NM_017304
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTATCCGGAATCAACCACGGGTCTCCAGCTCGACTCTCCCTGCGGCAGACAGGCTCCCCGGGATGA
TCTACAGTACTCGGTATGGGAGTCCAAAAGACAGCTCCAGTTTTACAGGAATCTGGGCAAATCTGGCCT
TCGGGTCTCCTGCCTGGGGCTTGAACATGGGTGACCTTCGGAGGCCAGATCACAGATGAGATGGCAGAG
CACCTAATGACCTTGGCCTATGACAATGGCATCAACCTGTTTCGATACGGCGGAGGTCTACGCAGCTGGCA
AGGCTGAAGTGGTATTAGGGAACATCATTAAGAAGAAGGGGTGGAGACGGTCCAGCCTTGTTCATCACCAC
CAAGATCTTCTGGGGCGAAAGGCAGAGACCGAGAGAGGCCTTCCCGGAAGCACATAATCGAAGGACTG
AAAGCTTCCCTGGAGAGGCTGCAGCTGGAGTACGTGGATGTGGTTTTTGCACACCGCCAGACCCCAACA
CACCCATGGAAGAGACTGTGCGGGCCATGACCCATGTCATCAACCAAGGGATGGCCATGTACTGGGGCAC
ATCAGCTGGAGCTCCATGGAGATCATGGAGGCCACTCGGTGGCTCGGCAGTTCAACCTGATCCCGCCC
ATCTGCGAGCAAGCGGAATACCACATGTTCCAGAGGGAGAAGGTAGAGGTCCAGCTGCCAGAGCTGTTCC
ACAAGATAGGAGTAGGTGCCATGACCTGGTCCCCTCTGGCCTGCGGCATTGTCTCAGGGAAGTATGACAG
TGGCATCCACCCACTCCAGAGCCTCTCTGAAGGGCTACCAGTGGTTGAAGGACAAGATCCTGAGTGAG
GAGGGTCCCGCCAGCAAGCCAAACTGAAGGAAGTGCAGGCCATTGCAGAGCGCCTAGGCTGCACCCCTCC
CACAGCTGGCCATAGCCTGGTGCCTGAGGAACGAGGGCGTCAGCTCCGTGCTCCTGGGTGCTTCCAATGC
AGAACAACCTTATGGAGAACATTGGAGCAATACAGGTCTTCCAAAATTGTCGTCCTCCATCGTCCACGAG
ATCGACAGCATTCTGGGCAATAAACCTACAGCAAAAAGGACTATAGATCC

ACGCGTACGCGGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)


Protein Sequence: >RR215276 representing NM_017304
 Red=Cloning site Green=Tags(s)

MYPESTTGSPARLSLRQTGSPGMIYSTRYGSPKRQLQFYRNLGKSGLRVSLGLGTWVTFGGQITDEMAE
 HLMTLAYDNGINLFDTAEVYAAGKAEVVLGNI IKKKGWRRSSLVITTKIFWGGKAETERGLSRKHIEGL
 KASLERLQLEYVDVVFANRPDPNTPMEETVRAMTHVINQGMAMYWGT SRWSSMEIMEAYSVARQFNLI PP
 ICEQAEYHMFQREKVEVQLPELFHKIGVGAMTWSPLACGIVSGKYDSGIPPYSRASLKG YQWLKDKILSE
 EGRRQQA KLKELQAIAERL GCTLPQLAI AWCLRNEGVSSVLL GASNAEQLMENIGAIQVLPKLS SIVHE
 IDSILGNKPYSKKDYRS


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_017304

ORF Size:	1101 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_017304.2 , NP_059000.1
RefSeq Size:	1590 bp
RefSeq ORF:	1104 bp
Locus ID:	29738
UniProt ID:	P62483
Cytogenetics:	5q36
MW:	41 kDa
Gene Summary:	voltage gated potassium channel beta subunit; modifies the inactivation kinetics of the channel [RGD, Feb 2006]