

Product datasheet for **RR215363**

Atg4c (NM_001107948) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Atg4c (NM_001107948) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Atg4c
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR215363 representing NM_001107948
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAGGCGTCAGGAACAGATGAAGTTGACAAGCTAAAACTAAATTTATATCTGCTTGAACAACATGA
AATATAGTTGGGTGTTGAAAACAAAGACATATTTAGTAGAACTCACCTGTTTTATTGCTTGGAAAATG
TTATCATTAAATATGAAGATGAAAGCAAGGTGTTGCCCGCAAGATCAGGATGTGCCATTGAAGATCAT
GTAATTGCAGGCAATGTGGAAGAGTTCGGTAAGGATTTCAATTTCTAGAATATGGCTGACCTACAGGGGAA
AATTTCTCAAATAGAAGCTTCAGCCTTACGACAGACTGTGGCTGGGGCTGTACTGCGGACCGGCCA
GATGCTGCTGGCTCAGGGACTCATACTACATTTCTCGGTAGAGCCTGGACCTGGCCGACGCACTGCAC
ATTGAGAGTTCGGACTCCGACTCATGGACCTCCAACACCATCCACAAGTTCAGTGCCTTTGAAGCCT
CCCTTTCAGGGGAGAGAGAACTCCGACCCAGCAGTCTCTCTGAAGGAAACGTCCGGGAAGCATCCCGA
CGACCATGCAGTACAGAGTGAGATTTACCACAGGCAGATCATCTTGGTTTGGAGACTCCCTGTGGCT
GTCTTTGGCTTACATCGCTGATAGAGTTCGGGAAGAAGTCTGGGAAGAAAGCTGGTACTGGTATGGAC
CAGCAGTGGTTGCTCACATTTAAAGAAAAGCAGTTGAAGAAGCAAGGCATCCGGATTTGCAAGGGCTAAC
AATTTATGTTGCCAGGACTGTACAGTTTACAATTCGATGTAATTGATAAGCAGACTGATTCTGTGACG
GCTGGAGATGCTCGTGACAAGGCTGTTATTATTTAGTTCCTGTTAGGCTTGGTGGAGAAAGAACCAACA
TTGACTATCTGGAGTTTGTGAAGGGTGTTTAAGCCTTGAGTATTGTGTTGGTATCATTGGCGGCAAAC
GAAACAGTCATACTACTTTGCTGGATTTCAAGATGACAGTTTGATTTACATGGATCCTCATTACTGCCAG
TCTTTGTAGATGTCAGCATAAAGGATTTCCCTCTCGAGACATCCACTGTCCATCTCCTAAGAAGATGT
CATTTGAAAAATGGATCCTAGTTGCACGATAGGATTTACTGTGAAAATGTCCAGGACTTTGAACGAGC
ATCTGAAGAAATCACTAAGATGTTGAAAATTTCTCTAAGGAGAAATATCCATTATTTACTTTTGTAAAT
GGCCATTCCAGAGACTTTGATTTACGTCCACTGCAGCCAGTGAGGAAGACCTTTTCTTAGAAGATGAAA
AGAAGAATTTCAAAGATTTAGCACAGAAGAGTTTGTCCCTCT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR215363 representing NM_001107948
Red=Cloning site Green=Tags(s)


MEASGTDEVDKLTKFISAWNNMKYSWVLKTKTYFSRNSPVLLLKGCYHFKEYEDESQVLPARGCAIEDH
 VIAGNVEEFRKDFISRIWLTYYEEFPQIEASALTTDCGWGCTLRGTQMLLAOGLILHFLGRAWTWPDALH
 IESSSDSWSNTIHKFTASFEASLSGERELRTPAVSLKETSQKHPDDHAVQSEIYHRQIISWFGDSPVA
 VFGLHRLIEFGKKSQKAGDWDYGPVVAHILRKAVEEARHPDLQGLTIYVAQDCTVYNSVDIKQTDSTV
 AGDARDKAVIILVPVRLGGERTNIDYLEFVKGVLSLEYCVGIIGGKPKQSYFFAGFQDDSLIYMDPHYCQ
 SFVDVSIKDFPLETFHCPSPKMSFRKMDPSCTIGFYCRNVQDFERASEEITKMLKISSKEKYPLFTFVN
 GHSRDFDTSTAASEEDLFLDEKKNFKRFSTEEFVLL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN:	NM_001107948
ORF Size:	1374 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001107948.1 , NP_001101418.1
RefSeq Size:	2331 bp
RefSeq ORF:	1377 bp
Locus ID:	313391
Cytogenetics:	5q33
MW:	52.1 kDa