

Product datasheet for **RR215364L3V**

Ngb (NM_033359) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ngb (NM_033359) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Ngb
Synonyms:	Glnh
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_033359
ORF Size:	453 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR215364).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_033359.3 , NP_203523.2
RefSeq Size:	1773 bp
RefSeq ORF:	456 bp
Locus ID:	85382
UniProt ID:	Q99JA8
Cytogenetics:	6q31
Gene Summary:	a porphyrin-containing protein; involved in oxygen supply of the neuron [RGD, Feb 2006]



[View online »](#)