

## Product datasheet for **RR215403**

### **Nfatc2 (NM\_001107805) Rat Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	Nfatc2 (NM_001107805) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Nfatc2
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide  
Sequence:

>RR215403 representing NM\_001107805  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCCCGGATCGCC

ATGGACGCCCCGGAGCCGAGTCCGACCCCGACGGCGGGGACGGCCCCGGCCACGAGCCGGGGCAGTT  
 CCCAAGACGAGCTGGACTTTCCATCCTCTTCGATTATGACTATCTGAACCTATCGAAGAAGAACCAT  
 CGCACATAAGGCCATCAGCTCACCTCCGGACTCGCATACCCGGATGATGTCCTGGACTATGGCCTCAAG  
 CCATGCAACCCTCTCGCCAGTCTCTCTGGCGAGCCCCCTGGCCGGTTCGGAGAGTCGGATAGTATAGGGT  
 TCCAGAACTTTCTGAGCCCGCCAAGCCAGCAGGGGCTTCGGGCCGAGCCCTCGGATCGAGACTCTCC  
 ATCCCACGAACTGATGCAGGCAGGGGGGCCCTCCGTGGGAGAGATGCCGGCCTGTCCCCGAGCAGCCG  
 GCCCTGGCTCTGGCCGGCATGGCCGCCAGCCGAGGTTACACTGCCCGTGCCGGCTTCGAGGGCTACC  
 GCGAGCCGCTTTGCTTGAGCCCCGCTAGCAGCGGCTCCTCTGCCAGTTCATTTCTGACACCTTCTCCCC  
 CTACACCTCGCCCTGCGTCTCACCAATAACGCCGGGCCGACGACCTGTGTCCCAGTTTCAAACATC  
 CCTGCTCATTATCCCCCAGAACCCTCTCCAATAATGTCACCTCGAACCAGCCTCGCCGAGGACAGCTGCC  
 TGGGCCGACACTCGCCCGTGCCCGTCCGGCATCCCGCTCCTCCTCACCCGGTGCCAAGCGGAGGCATTC  
 GTGCGCAGAGGCCCTTGTTGCCCTCTGCCCGCAGCCTCACCCAGCGCTCCCGGAGCCCTCGCCACAG  
 CCCTCGCTCAGTGGCACCCACAGGACGACGGCATCCCGCTGGGTACCCCCACGGCCGGCTCTGCTG  
 TTCTCATGGATGCCCTCAACAGCCTGGCCACCGACTCGCCCTGCGGGATCCCCTAAGATATGGAAGAC  
 CAGTCTGATCCAACCTCCCGTGTCCACTGCTCCATCCAAGGCAGGCCTGGCCGCCACATCTACCTGT  
 GTGGAGTTCCTGGGGCCGTGTGAGCAGGAGGAGAGGAATTCGCTCCAGAGTCCATCCTGTGGTAC  
 CACCCACTGGCCCAAGCAATTGGTGGCCGCAATCCCATCTGCAGCATCCCCGTGATGCATCCCTCCC  
 ACCACTAGAGTGGCCACTCTCAACCAAGTCCGGCTCCTATGAGCTTCGATTGAGGTCCAGCCCAAGCCC  
 CATCACCGAGCCCACTATGAGACAGAGGCGAGCCTGGCGCTGTCAAAGCCCCAACAGGAGGACACCCTG  
 TGGTGCAGCTCCACGGCTACATGGAGAACAAGCCTCTGGGGTTCAGATCTTCATCGGGACGGCTGATGA  
 ACGGATCCTTAAGCCACACGCTTCTACCAAGTACACAGGATCACTGGGAAAACAGTGACAACCACTAGC  
 TATGAGAAGATTGTGGGCAACACCAAGTCTGGAGATCCCCCTGGAGCCAAAGAACAACATGCGGGCCA  
 CCATCGACTGCGCAGGCATCCTGAAGCTCCGAAACGCTGACATCGAGCTGCGGAAGGGTGAACAGACAT  
 AGGCAGAAAGAACACGCGTGTGCGCCTGTGTTCCGTGTGCACATCCCGGAGCCAGTGGACGCATCGTC  
 TCCTGCAGGCTGCATCAACCCCATCGAGTGTCTCAGCGCTCTGCCACGAGCTGCCATGGTGGAGA  
 GACAAGACATGGACAGCTGCCTGGTCTACGGCGCCAGCAGATGATCCTCACGGCCAGAACTTCACAGC  
 TGAGTCCAAGGTTGTGTTTCATGGAGAAGACTACAGATGGGCAGCAGATTTGGGAGATGGAAGCTACAGT  
 GATAAAGACAAGAGCCAGCCTAACATGCTTTTTGTTGAGATCCCCGAGTATCGGAACAAGCACATCCGGG  
 TGCCCGTGAAGTGAACCTTCTACGTCATCAACGGAAAGAGAAAACGAAGTCAAGCCACAGCACTTTACTTA  
 CCACCCAGTCCCCGCCATCAAGACAGAGCCAGCGATGAGTATGAACCGTCTTTGATCTGCAGCCCCGCC  
 CACGGAGGCTGGGGAGCCAGCCATATTACCCACAGCACCAATGCTGGCCGAGTCCCCCTCCTGTCTTG  
 TGGCCACCATGGCCCCTTGCCAACAGTTCGGCTCAGGGCTCTCATCCCCTGATGCTCGTACCAGCAGCA  
 GAGCCCCGACGCCCTCTACCAGAGGAGCAAGAGCCTGAGTCCAGGCCTGCTGGGCTACCAGCAGCCG  
 TCCCTCTTGGCAGCACCCCTGGGTCTGGCTGATGCCACCGTTCGTGCTGGTGCATGCTGGTTCGAGG  
 GGCAGGGGACGGGACGGGCTCCACCTTCCACACACCTCCTCGGCCAGCCAGCAGGCCTCACCCGTGAT  
 CCACTACTCACCCACCAACCAGCAGTTCGTGGTGGGGTCAACAGGAGTTCAGCATATCATGACTGT  
 GAGAATTCGGCCCCAGCTCTGTGAGGCTGGCCTGCCTCCCGTCAACCAAGGTGAGAGGCTGAGCCAG  
 GGTCTACCCACGGTCACTCAACAACAGACTGCCCGAGCCAAAGAGCTGCCAAAACGGACCCCGGT  
 CAGTGACCAGAAGGAAGCCCTGCCAGTGGAGTACAGTCAAACAGGAGCAGAACCTGGACCAGACTTAC  
 CTGGATGACGTTAATGAAATCATCAGGAAGGAGTTTTCCGGACCTCCCTCCGAAATCAGACC

AGCGGACCGACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC  
 TGGATTACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RR215403 representing NM\_001107805  
 Red=Cloning site Green=Tags(s)


MDAPEPQSDPDGGDGPHEPGGSSQDELDFILFDYDYLNPIDEEPIAHKAISSPSGLAYPDDVLDYGLK  
 PCNPLASLSGEPPGRFGESDSIGFQNFLSPAKPAGASGSPRIETPSHELMQAGGALRGRDAGLSPEQP  
 ALALAGMAASPRFTLPVPGFEGYREPLCLSPASSGSSASFISDTFSPYTPSPCVSPNAGPDDLCPQFQNI  
 PAHYSPRTSPIMSPRTSLAEDSCLGRHSPVPRPASRSSSPGAKRRHSCEALVAPLPAASQRSRSPSPQ  
 PSHVAPQDDGIPAGYPPTAGSAVLM DALNSLATDSPCGIPSKIWKTS PDTPVSTAPSKAGLARHIYPA  
 VEF L G P C E Q E E R R N S A P E S I L L V P P T W P K Q L V P A I P I C S I P V T A S L P P L E W P L S N Q S G S Y E L R I E V Q P K P  
 H H R A H Y E T E G S R G A V K A P T G G H P V V Q L H G Y M E N K P L G L Q I F I G T A D E R I L K P H A F Y Q V H R I T G K T V T T T S  
 Y E K I V G N T K V L E I P L E P K N M R A T I D C A G I L K L R N A D I E L R K G E T D I G R K N T R V R L V F R V H I P E P S G R I V  
 S L Q A A S N P I E C S Q R S A H E L P M V E R Q D M D S C L V Y G G Q Q M I L T G Q N F T A E S K V V F M E K T T D G Q Q I W E M E A T V  
 D K D K S Q P N M L F V E I P E Y R N K H I R V P V K V N F Y V I N G K R K R S Q P Q H F T Y H P V P A I K T E P S D E Y E P S L I C S P A  
 H G L G S Q P Y P Q H P M L A E S P S C L V A T M A P C Q Q F R S G L S S P D A R Y Q Q Q S P A A A L Y Q R S K S L S P G L L G Y Q Q P  
 S L L A A P L G L A D A H R S V L V H A G S Q G Q G Q G S T L P H T S S A S Q Q A S P V I H Y S P T N Q Q L R G G G H Q E F Q H I M Y C  
 E N F G P S S V R P G L P P V N Q Q R L S P G S Y P T V I Q Q T A P S Q R A A K N G P P V S D Q K E A L P S G V T V K Q E Q N L D Q T Y  
 L D D V N E I I R K E F S G P P S R N Q T


SGPTRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:**

SgfI-RsrII

**Cloning Scheme:**



**Plasmid Map:**


**ACCN:** NM\_001107805

**ORF Size:** 2793 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001107805.1](#), [NP\\_001101275.1](#)

**RefSeq Size:** 3701 bp

**RefSeq ORF:** 2796 bp

**Locus ID:** 311658

**Cytogenetics:** 3q42

**MW:** 100.3 kDa