

Product datasheet for **RR215430**

Npy1r (NM_001113357) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Npy1r (NM_001113357) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Npy1r
Synonyms: NPY-1
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR215430 representing NM_001113357
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGAATCAACTCTGTTCTCCAGGGTTGAAAATTACTCAGTTCCTATAATGTCTCAGAGAATTCTCCAT
TTCTGGCCTTTGAGAATGATGACTGCCACCTGCCCTTGGCTGTGATTCACCTTAGCTCTTGCTTATGG
GGCTGTGATTATCTTGGGGTCTCTGAAACCTGGCATTGATCATAATCATCCTAAAACAGAAAGAGATG
AGGAATGTACCAACATTCTGATCGTGAACCTCTCCTTCTCAGACTTGCTGGTTCGAGCATGTGTCTCC
CGTTACCTTTGTGTACACTGATGGACCACTGGGTCTTCGGGGAGACCATGTGCAAACTGAATCCTTT
TGTGCAATGCGTCTCCATTACAGTATCCATTTCTCTGTTTCTCATCGCTGTGGAACGTATCAGCTA
ATCATCAACCAAGAGGGTGGAGACCAACAATAGACATGCTTACATAGGCATTACTGTCATTTGGGTAC
TGGCGGTGGCTTCTTCTGCCCCGTCGTGATCTATCAAATTCGACTGATGAGCCCTTCCAAAATGTATC
ACTTGGCGGTCAAGGACAAGTATGTATGTTTGACAAATCCCGTCAGACTCTCACAGGCTGTCTTAC
ACGACTCTTCTTCTGGTGTGCAGTATTTGGCCACTCTGTTTCATATTCATATGCTACTTCAAGATAT
ACGTTTCGCTTAAAAGGAGAAACAACATGATGGACAAAATCAGGGACAGTAAGTACAGTCCAGTGAGAC
CAAACGAATCAACGTCATGCTGCTCTCCATCGTGGTTCGCTTTCGCGGTCTGCTGGCTGCCCTTACCATC
TTCAACACTGTGTTTCGACTGGAATCACCAGATCATTGCCACCTGCAACCAATCTGCTGTTCTGCTCT
GCCACCTCACGGCCATGATCTCCACCTGCGTCAACCCCATCTTTTATGGATTCTGAACAAAAATTTCCA
GAGAGACTTGAGTTCTTCTTCACTTTTGTGACTTCCGGTCTCGAGACGACGACTACGAGACTATAGCC
ATGTCTACCATGCATACGGACGTGTCCAAGACTCTTTGAAGCAGGCAAGCCCGTTGCATTTAAAAAA
TCAGTATGAATGACAATGAAAAATC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR215430 representing NM_001113357
Red=Cloning site Green=Tags(s)


MNSTLFSRVENYSVHYNVSENSPFLAFENDDCHLPLAVIFTLALAYGAVIILGVSGNLALIIIIILKQKEM
 RNVNTNILIVNLSFSDLLVAVMCLPFTFVYTLMDHWVFGETMCKLNPVQCVSITVSIFSLVLI AVERHQL
 IINPRGWRPNRRHAYIGITVIWVLAVASSLPFVIYQILTDEPFQNVSLAAFKDKYVCFDKFPSSDHRLSY
 TLLLLVLQYFGPLCFIFICYFKIYVRLKRRNNMMDKIRD SKYRSSETKRINVMLLSIVVAFVCWLP LTI
 FNTVFDWNHQIIATCNHLLFLLCHLTAMISTCVNPIFYGFLNKNFQRDLQFFNFCDFRSRDDDYETIA
 MSTMHTDVSKTSLKQASPVAFKKISMNDNEKI

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001113357

ORF Size:	1146 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001113357.1 , NP_001106828.1
RefSeq Size:	2618 bp
RefSeq ORF:	1149 bp
Locus ID:	29358
Cytogenetics:	16p14
MW:	44.1 kDa
Gene Summary:	G protein-coupled receptor for neuropeptide Y and peptide YY [RGD, Feb 2006]