

Product datasheet for **RR215445L3V**

Ube2d3 (NM_031237) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ube2d3 (NM_031237) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Ube2d3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_031237
ORF Size:	441 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR215445).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_031237.1 , NP_112516.1
RefSeq Size:	1531 bp
RefSeq ORF:	444 bp
Locus ID:	81920
UniProt ID:	P61078
Cytogenetics:	2q43
Gene Summary:	enzyme that has an important role in ubiquitin mediated proteolysis by conjugating ubiquitin to protein substrates [RGD, Feb 2006]



[View online »](#)