

Product datasheet for **RR215452L3V**

Fgf5 (NM_022211) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Fgf5 (NM_022211) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Fgf5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_022211
ORF Size:	798 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR215452).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_022211.1 , NP_071547.1
RefSeq Size:	801 bp
RefSeq ORF:	801 bp
Locus ID:	60662
UniProt ID:	P48807
Cytogenetics:	14p22
Gene Summary:	member of the fibroblast growth factor (FGF) family; implicated in activities including embryonic development, cell growth, morphogenesis, tissue growth, and tumor growth and invasion [RGD, Feb 2006]



[View online »](#)