

Product datasheet for **RR215472**

Med15 (NM_001108325) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Med15 (NM_001108325) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Med15
Synonyms:	Pcqap
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)


ORF Nucleotide
Sequence:

>RR215472 representing NM_001108325
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGACGTTTCGGGGCAGGAGACCGACTGGCGTAGCGCCGCTTTCGGCAGAAGCTGGTCAGCCAAATTG
AGGATGCCATGAGGAAAGCTGGTGTGGCCACAGTAAATCTAGCAAGGATATGGAGAGTCATGTGTTCTT
GAAGGCCAAGACCCGGGATGAATATCTTCCCTTGTGGCCGACTCATTATCCATTTCCGAGATATTCAT
AACAGAAGATCCCAAGCTTCTGTGACTGACCCATGAATGCACTGCAGAGCCTTACTGGTGGACCCACCC
CAGGAGCAGCTGGGATTGGCATGCCTCTCGGGGCCAGGACAGTCCCTGGTGGGATGGTGGCCTTGG
TGCGATGGGACAACCACTGCCCTCTCCGGGCAGCCACCCCTGGAACCTCTGGAATGGCCCTCATGGC
ATGGCTGTGGTGTCTACAGCACTCCACAGACTCAGCTGCAGTCCAGCAAGTGGCATTGCAGCAGCAGC
AACAGCAACAACAGCAGCAGTTCAGCAACAGCAGGCAGCATTGCAGCAGCAGCAACAGCAGCAGCAGCA
ACAGCAGCAGCAACAGCAGTTCAGGCACAACAGAATGCCATGCAGCAGCAGTTCGAAGCAGTAGTGCAG
CAGCAGCAACTACAGCAGCAGCAGCAGCAACAGCACCTGATTAAGTTGCATCACCAAGCCAGCAACAAC
AGATGCAACAGCAGCAGCAACTGCAGAGGATGGCAGATTGCAGTGCAGCAACAGCAACAACAGCAGCA
ACAGCAGCAGGCTTTCAGGCCCAGCCACCAATGCAACAGCCACCAATGCAGCAGCCACAGCCTCTCTCT
TCTCAGGCCCTACCCAGCAGCTGTACAGATGCATCATCCACAGCATCACCCAGCCACCCTCAGGCTC
AGCAGTCCCCCGTTGCTCAGAACCAACCACCAACTCCAACCACAGTACAGAACCCAGCCTTTGGTGTG
ACAAGCACAAGCCCTGTCTGGACCAATGCTCTATGCTGCCAACAGCCACTGAAATTTGTCGGTCTCCG
ATGGTGGTCCAGCAGCCGAGGTGCAGCCCAGGTGCAGCAGGTGCAGCCCAGGTGCAGCCACAGGCAG
ACGTGCAGACAGCAGCAGTCTGCCAGATGGTAGCTCCCGGCTTCAGATGATTGCTGAAGCCTTGGCCCA
AGGCGGGATGCACGTAAGAGCCCGTTCCCGCCACCTCCACCATGTCTGCTGGCCCGTCAAGCTCCATC
TCTTTGGCGGCAGCCACGGCACAGGTGAGCCAAAGCAGCCTCACCATGCTGTCTCACCGTCCACAG
GCCAGCAGGTGCAGACCCACAGTGCATGCCCTTCCCCACAGCCATCCCCACAACCTGGCTCACAGCC
CAACTCTAATGTGAGCTCCGGCCCTGCCCATCCCCAGCAGTTCCTGCCTAGCCCTCACACAGCCT
TCCAGAGCCAGTGACAGCAGCAGCCACAGAACTTCAGCGTCCCTTCCCCGGACCTTTAACACCC
CTGTGAACCCAGCTCTGTGATGAGCCCGCTGGCTCCAGTCCAGTCCAGTCCAGTCCAGTCCAGTCCAGT
GCTGAAGCAACTGTCCAAGTACATCGAGCCCTTGCCTGATGATCAACAAGATCGACAAGAACGAAGAC
AGGAAAAAGGACTTAAGTAAGATGAAGAGCTTGTGGACATCCTCACTGACCCCTTAAGCGGTGTCCC
TGAAGACCCTGCAAAAGTGTGAGATTGCCCTGGAGAAGCTCAAGAATGACATGGCAGTGCCACACCCCC
ACCACCCCGATTCTTCAACCAAAACAGCAGGACCTGTGCCAACCACTCCTAGATGCAGTCCAGGCAAC
ATCCGTTACCTGTCTTCAACATTCCCTGTATCGCACATTCGTGCCAGCCATGATGGCCATCCATGGTC
CGCCTATCGTGTCCCAAGTGGTGTCTCCCGGAAGCGCAGGTTTGGAGGAGATGAGCGGCAGAGCATACC
CAATGTGTGCAAGGTGAGGTGGCAAGGCTGGATCCCAAGTTCCTGGTGAACCTGGACCTTCTCACTGC
AGCAACAACGGTACTGTCCACCTGATCTGCAAGCTGGATGACAAGGACCTCCCTAGTGTGCCACCACTGG
AGCTCAGTGTACCTGTGACTACCCTGCCAGAGCCCGATGTGGTTCGACCGTCAAGTGGCAATATGATGC
CAACCCCTTCTGCAGTCAAGTGCACCGGTGCATGACCTCCAGGCTGCTGCAGTCCCTGACAAGCACTCA
GTCACAGCCCTGCTCAACACCTGGGCCAGAGCATCCACAGGCTGCTTCTCAGTGGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Plasmid Map:


ACCN: NM_001108325

ORF Size: 2370 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001108325.1](#), [NP_001101795.1](#)

RefSeq Size: 3243 bp

RefSeq ORF: 2373 bp

Locus ID: 360743

Cytogenetics: 11q23

MW: 86.8 kDa