

## Product datasheet for **RR215474L3V**

### Siah1 (NM\_080905) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Siah1 (NM_080905) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Siah1
Synonyms:	SIAH; Siah1a
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_080905
ORF Size:	846 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR215474).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_080905.2</a> , <a href="#">NP_543181.2</a>
RefSeq Size:	2105 bp
RefSeq ORF:	849 bp
Locus ID:	140941
UniProt ID:	<a href="#">Q920M9</a>
Cytogenetics:	19p11
Gene Summary:	synaptophysin binding protein and function as E3 ubiquitin-protein ligases to regulate the ubiquitination and degradation of synaptophysin [RGD, Feb 2006]



[View online »](#)