

Product datasheet for **RR215477L4V**

Slc2a4 (NM_012751) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Slc2a4 (NM_012751) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Slc2a4
Synonyms:	Glut4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_012751
ORF Size:	1527 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR215477).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_012751.1 , NP_036883.1
RefSeq Size:	2506 bp
RefSeq ORF:	1530 bp
Locus ID:	25139
UniProt ID:	P19357
Cytogenetics:	10q24
Gene Summary:	transports glucose; may play a role in insulin response [RGD, Feb 2006]



[View online »](#)