

Product datasheet for **RR215486**

Slc6a11 (NM_024372) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Slc6a11 (NM_024372) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Slc6a11
Synonyms:	Gabt4; Gat3
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>RR215486 representing NM_024372
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGACTGCGGAGCAAGCGCTGCCCTGGGCAACGGAAGCGGCCGAGGAGCGCGAGGGTCCGAGGCGC
 TGGGCGGCGGCGGGGCGCGCGGGACGCGAGGCGCGACAAGGCGGTCCACGAGCGCGGTCA
 CTGGAACAACAAGGTGGAGTTCGTGTTGAGCGTAGCGGGAGAGATCATCGGTCTGGGCAACGTGTGGCGC
 TTCCCCTACCTGTGCTACAAGAACGGCGGAGGGCATTCTGATTCTTACGTGGTGTTTTTTTCATCTGCT
 GTGGAATCCCCGCTTCTTCTGAAACGGCTCTGGGCGAGTTCACGAGCGAGGGCGGCATCACGTGCTG
 GAGGAGAGTCTGCTTTATTTGAAGGCATCGGCTATGCAACACAGGTGATCGAGGCGCATCTCAATGTC
 TACTACATCATCATCTGGCGTGGGCCATCTTCTACTTAAGCAACTGCTTACCACCGAGCTCCCCTGGG
 CCACCTGTGGGCATGAGTGAACACAGAGAAATGTGTGGAGTTCAGAAGCTGAAGTTCAGCAACTACAG
 TCATGTGTCCCTGCAGAACGCAACTCCCCTGATGGAGTTCGGAACGCCGGTCTTGCTATATCT
 GATGGCATTGAACACATCGGGAACCTCCGATGGGAGCTGGCACTGTGTCTCTGGCGGCTTGACCATCT
 GCTACTTCTGCATCTGGAAGGTACGAAGTCAACTGAAAAGTCTGTATGTCACCTGCAACCTTCCCCTA
 CATCATGCTGCTGATCCTCCTGATCCGAGGGGTACGTTGCCGGTGCCTCGGAAGGCATCAAGTCTAC
 CTGTACCTGACCTCTCCCGGCTCTGATCCACAGGTGTGGGTGGATGCTGGGACGCAGATCTTTTTCT
 CCTATGCCATCTGCCTGGGCTGCCTGACCGCTCTGGGGAGTTACAACAATAACAACAAGTCTACAG
 GGACTGTATTATGCTCTGCTGTGAACAGTGGCACCAGTTCGTGGCTGGGTTTGTATCTTCTCAGTC
 CTGGGCTTCATGGCGTACGAGCAGGGCGTGCCTATTGCTGAGGTGGCAGAATCAGGTCTGGACTGGCTT
 TCATCGCCTACCCCAAGGCTGTCACTATGATGCCCTGTCCCATTTGTGGCCACCCTGTTCTTCATGAT
 GCTCATCTTCTGGCCCTGGACAGTCAAGTTTGTGTGTGGAGAGCCTTGTGACAGCCGTGGTTGACATG
 TACCCCAAGTCTTCCGGCGGGCTACCGGCGAGAAGTCTCATCTGGCCCTGTCCATTGTCTCTTATT
 TCCTAGGCCTGGTGTGCTGACAGAGGGAGGCATGTACATTTTCCAGCTTTTGTACTCATACGCCCCAG
 TGGCATGTGCTTCTCGTGGCCATCTTGTAGTGTCTGCATCGGCTGGGTGTATGGAAGTAACAGG
 TTCTATGACAATATTGAGGACATGATTGGATACCGGCCACTGCTACTCATCAAGTGGTGTGAAAGTTG
 TGACCCCTGGGATCTGTGCGGGCATCTTCATCTTCTTCTGGTCAAGTACAAGCCGCTCAAGTACAACA
 TGTGTACACATATCTGCTTGGGCTACGGCATTGGCTGGCTCATGGCTGTCTCCATGCTGTGCATC
 CCGCTCTGGATCTTCATCAAGCTGTGGAAGACAGAGGGCACCCCTGCCCGAGAAATTACAGAAAGTTGACAG
 TCCCAGCGCTGATCTGAAAATGAGGGCAAGCTTGGGGCCAGCCACCGATGGTACCCTTAATGACTG
 TGAGGCCAAGGTCAAAGGCGACGTACCATCTCTGCCATCACAGAGAAGGAGACGCACTT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:


>RR215486 representing NM_024372
 Red=Cloning site Green=Tags(s)

MTAEQALPLGNGKAAEEARGSEALGGGGGAAGTREARDKAVHERGHWNKVEFVLSVAGEIIGLGNVWR
 FPYLCYKNGGGAFLIPYVVFICCGIPVFFLETALGQFTSEGGITCWRRCVPLFEGIGYATQVIEAHLNV
 YYIIILAWAIFYLSNCFTELWPATCGHEWNTKCVFQKLNFSNYSHVSLQNATSPVMEFWERRVLAIS
 DGIEHIGNLRWELALCLLAAWTICYFCIWKGTSTGKVVYVTAATFPYIMLLILLIRGVTLPGASEGIKFY
 LYPDLSRLSDPQVVDAGTQIFFSYAICLGLTALGSYNNYNNCYRDCIMLCCLNSGTSFVAGFAIFSV
 LGFMAYEQVPIAEVAESGPGFLAIYPKAVTMMPLSPLWATLFFMMLIFLGLDSQFVCVESLTVAVVDM
 YPKVFRRGYRRELLILALSIVSYFLGLVMLTEGMYIFQLFDSYAASGMCLLFVAIFECVICGWVYGSNR
 FYDNIEDMIGYRPLSLIKWCWKVVTGICAGIFIFLVKYKPLKYNVYTYPAWGYGIGWLMALSSMLCI
 PLWIFIKLWKTEGTLPEKLQKLTVPADLKMGRKLGASPRMVTVNDCEAKVKGDGTISAITEKETHF


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_024372

ORF Size: 1881 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_024372.2 , NP_077348.1
RefSeq Size:	1935 bp
RefSeq ORF:	1884 bp
Locus ID:	79213
UniProt ID:	P31647
Cytogenetics:	4q42
MW:	69.9 kDa
Gene Summary:	catalyzes sodium and chloride dependent transport of gamma-aminobutyric acid (GABA) to terminate GABA mediated synaptic activity [RGD, Feb 2006]