

Product datasheet for **RR215490L4V**

Rpl32 (NM_013226) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Rpl32 (NM_013226) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Rpl32
Synonyms:	MGC72905
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_013226
ORF Size:	405 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR215490).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_013226.2 , NP_037358.1
RefSeq Size:	564 bp
RefSeq ORF:	408 bp
Locus ID:	28298
UniProt ID:	P62912
Cytogenetics:	4q42
Gene Summary:	ribosomal protein subunit [RGD, Feb 2006]



[View online »](#)