

Product datasheet for **RR215494**

Sftpc (NM_017342) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Sftpc (NM_017342) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Sftpc
Synonyms: SP-C
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR215494 representing NM_017342
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGACATGGGTAGCAAAGAGGTACTGATGGAGAGCCCACGGATTACTCGACAGGTCCCAGGAGCCAGT
TTCGATTCCCTGCTGCCCGTGCATCTCAAACGCCTTCTCATCGTGGTTGTGGTGTAGTCTTGTGCGT
CGTGGTATTGTAGGGGCTCTGCTCATGGGCCTTCACATGAGTCAGAAACATACTGAGATGGTCCTTGAG
ATGAGCATCGGAGGAGCACCGGAAACTCAGAAACGCCTCGCCCTGAGTGAACACACAGACACCATCGCTA
CCTTCTCTATTGGCTCCACTGGCATCGTTCTATATGACTACCAGCGGCTCCTGACCGCTATAAGCCAGC
TCCAGGAACCTACTGCTACATCATGAAGATGGCTCCAGAGAGCATCCCTAGTCTTGAGGCTTTGGCTAGA
AAATTCAGAAGCTCCAGGCCAAGTCTCCACACCCACCTCTAAGTTGGGCCAGGAAGAAGGGCACAGTG
CTGGATCCGATTCTGATTCTTCTGGGCGAGACCTGGCTTCTGGGCCTTGCTGTGAGCACCTTGTGTGG
AGAGCTACCACTGTACTATATC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR215494 representing NM_017342
Red=Cloning site Green=Tags(s)

MDMGSKEVLMESPPDYSTGPRSQFRIPCCPVHLKRLLI VVVVVVIVVVVIVGALLMGLHMSQKHTEMLVLE
MSIGGAPETQKRLALSEHTDTIATFSIGSTGIVLYDYQRLLTAYKPAPGTYYCYIMKMAPESIPSLEALAR
KFKNFQAKSSTPTSKLGQEEGHSAGSDSDSSGRDLAFLGLAVSTLCGELPLYYI

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_017342

ORF Size: 582 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_017342.2](#), [NP_059038.1](#)

RefSeq Size: 782 bp

RefSeq ORF: 585 bp

Locus ID: 50683

UniProt ID: [P11685](#)

Cytogenetics: 15p11

MW: 21 kDa

Gene Summary: promotes reduction of surface tension in the lung [RGD, Feb 2006]