

Product datasheet for RR215523L4

Ets1 (NM_012555) Rat Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ets1 (NM_012555) Rat Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Ets1
Synonyms:	Ets-1; Etsonc; Tpl1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR215523).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.



[View online »](#)

Locus ID: 24356
UniProt ID: [P41156](#)
Cytogenetics: 8q21
Gene Summary: activates transcription of cyclin-dependent kinase inhibitor p21WAF1/Cip1; plays a role in inhibition of apoptosis [RGD, Feb 2006]